

attys-comm

Generated by Doxygen 1.8.13

## Contents

<b>1</b>	<b>AttysCOMM</b>	<b>1</b>
<b>2</b>	<b>Hierarchical Index</b>	<b>5</b>
2.1	Class Hierarchy . . . . .	5
<b>3</b>	<b>Class Index</b>	<b>5</b>
3.1	Class List . . . . .	5
<b>4</b>	<b>Class Documentation</b>	<b>5</b>
4.1	AttysComm Class Reference . . . . .	5
4.1.1	Detailed Description . . . . .	6
4.1.2	Constructor & Destructor Documentation . . . . .	6
4.1.3	Member Function Documentation . . . . .	6
4.2	AttysCommBase Class Reference . . . . .	7
4.2.1	Detailed Description . . . . .	9
4.2.2	Member Function Documentation . . . . .	9
4.2.3	Member Data Documentation . . . . .	13
4.3	AttysCommListener Struct Reference . . . . .	20
4.3.1	Detailed Description . . . . .	20
4.4	AttysCommMessage Struct Reference . . . . .	20
4.4.1	Detailed Description . . . . .	20
	<b>Index</b>	<b>21</b>

## 1 AttysCOMM

The C++ & Python API for the Attys bluetooth data acquisition board: <http://www.attys.tech>

The library is cross platform: It's for Linux, Windows and Mac.

## Installation instructions

### Linux

```
cmake .  
make  
sudo make install
```

This will generate: a dynamic library libattyscomm\*.so and a static one called libattyscomm\_static.a.

### Windows:

Under windows only the static library is generated which should be used for your code development.

```
cmake -G "Visual Studio 15 2017 Win64" .
```

and then start Visual C++ and compile it.

### MacOS

For pure commandline install:

```
cmake .  
make  
make install
```

This will generate: a dynamic library libattyscomm.\*dylib and a static one called libattyscomm\_static.a.

If you want to debug/develop the library in Xcode:

```
cmake -G Xcode
```

## Usage

A small test program is in the `examples` directory which scans for an Attys and then prints the incoming data to stdout. Type `cmake .`, `make` and then `./attystest` to run it.

Here is a step by guide how to code it:

#### 1. scan for Attys

```
int ret = attysScan.scan();
```

#### 2. Check the number of Attys detected

```
attysScan.getNAttysDevices()
```

#### 3. If devices have been detected you can get them via `getAttysComm(0, 1, 2, etc)`.

#### 4. Set the parameters, for example:

```
attysScan.getAttysComm(0)->setAdc_samplingrate_index(AttysComm::ADC_RATE_250HZ);
```

#### 5. Register a callback (optional)

```
attysCallback = new AttysCallback(this);  
attysScan.getAttysComm(0)->registerCallback(attysCallback);
```

You need to overload the virtual function of the callback in your program.

#### 6. Start data acquisition

```
attysScan.getAttysComm(0)->start();
```

#### 7. Check if ringbuffer contains data and wait till true

```
attysScan.getAttysComm(n)->hasSampleAvilabale();
```

#### 8. Get samples from buffer

```
float* values = attysScan.getAttysComm(n)->getSampleFromBuffer();
```

#### 9. go back to 7)

#### 10. Ending the program:

```
attysScan.getAttysComm(n)->quit();
```

### Python (SWIG)

This library is fast and multi threaded. It performs the data acquisition in the background while python can then do the postprocessing.

Pre-compiled packages for Linux, Windows and MacOS are available.

#### Linux

##### Python package (pip):

Make sure you have the bluetooth development libraries:

```
sudo apt-get install libbluetooth-dev
```

and then install with:

```
pip3 install pyattyscomm
```

##### From source

You need to have swig-3.x installed. Then run:

```
cmake .  
make  
make install  
./setup.py install
```

and then you can load the module `pyattyscomm` system-wide!

## Windows

### Python package (pip):

In the python console type:

```
pip install pyattyscomm
```

### From source

Install `swig` and re-run the C++ installation. Make sure to select "Release" in Visual Studio as python is usually not installed with its debug libraries. After compilation you get:

- `Release\_pyattyscomm.exp`
- `Release\_pyattyscomm.pyd`
- `pyattyscomm.py`

Install them with:

```
python setup.py install
```

## MacOS

### Python package (pip):

```
pip3 install pyattyscomm
```

### From source

You need to have `swig-3.x` installed for homebrew. Then run:

```
cmake .
make
make install
./setup.py install
```

and then you can load the module `pyattyscomm` system-wide!

### How to use

The python API is identical to the C++ one. All the definitions are in [AttysComm.h](#) and `AttysScan.h`.

Here is an example:

```
# load the module
import pyattyscomm

# Gets the AttysScan class which scans for Attys via bluetooth
s = pyattyscomm.AttysScan()

# Scan for Attys
s.scan()

# get the 1st Attys
c = s.getAttysComm(0)

# if an attys has been found c points to it. Otherwise it's None.

# Start data acquisition in the background
c.start()

# Now we just read data at our convenience in a loop or timer or thread
while (not c.hasSampleAvilabale()):
    # do something else or nothing
    a = a + 1
    # getting a sample
    sample = c.getSampleFromBuffer()

    # do something with the sample
    print(sample)
```

## Demo programs

There are demo programs which show you how to read/plot data with pyattyscomm:

```
readdata_demo.py
realtime_plot_demo.py
realtime_two_channel_plot.py
```

Enjoy!

<http://www.attys.tech>

## 2 Hierarchical Index

### 2.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

<b>AttysCommBase</b>	<b>7</b>
<b>AttysComm</b>	<b>5</b>
<b>AttysCommListener</b>	<b>20</b>
<b>AttysCommMessage</b>	<b>20</b>

## 3 Class Index

### 3.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:


<b>AttysComm</b>	<b>5</b>
<b>AttysCommBase</b>	<b>7</b>
<b>AttysCommListener</b>	<b>20</b>
<b>AttysCommMessage</b>	<b>20</b>

## 4 Class Documentation

### 4.1 AttysComm Class Reference

```
#include <AttysComm.h>
```

Inheritance diagram for AttysComm:



## Public Member Functions

- [AttysComm](#) (void \*\_btAddr=NULL, int \_btAddrLen=0)
- virtual void [connect](#) ()
- virtual void [closeSocket](#) ()
- virtual void [start](#) ()
- unsigned char \* [getBluetoothBinaryAddress](#) ()
- void [getBluetoothAddressString](#) (char \*s)

## Additional Inherited Members

### 4.1.1 Detailed Description

[AttysComm](#) contains all the necessary comms to talk to the Attys on both Linux and Windows.

1) Instantiate the class AttyScan and do a scan It finds all paired Attys and creates separate [AttysComm](#) classes 2) These classes are in in the array attysComm in AttyScan and the number of them in nAttysDevices. 3) All attys↔Comm are Threads so just start the data acquisition with [start\(\)](#), for example attysComm[0]->[start\(\)](#) for the 1st Attys 4) Get the data either via the RingBuffer functions or register a callback to get the data as it arrives. [AttysComm](#) class which contains the device specific definitions and implements the abstract classes of [AttysCommBase](#). See [AttysCommBase](#) for the definitions there.

### 4.1.2 Constructor & Destructor Documentation

#### 4.1.2.1 AttysComm()

```
AttysComm::AttysComm (
    void *_btAddr = NULL,
    int _btAddrLen = 0 ) [inline]
```

Constructor: Win/Linux: takes the bluetooth device structure and its length as an argument. For Mac: just a pointer to the device.

### 4.1.3 Member Function Documentation

#### 4.1.3.1 closeSocket()

```
virtual void AttysComm::closeSocket ( ) [virtual]
```

closes socket safely

Implements [AttysCommBase](#).

#### 4.1.3.2 connect()

```
virtual void AttysComm::connect ( ) [virtual]
```

connects to the Attys by opening the socket throws exception if it fails.

Implements [AttysCommBase](#).

#### 4.1.3.3 getBluetoothAdressString()

```
void AttysComm::getBluetoothAdressString (
    char * s ) [virtual]
```

returns the Mac address as a string

Implements [AttysCommBase](#).

#### 4.1.3.4 getBluetoothBinaryAdress()

```
unsigned char* AttysComm::getBluetoothBinaryAdress ( ) [virtual]
```

returns an array of 14 bytes of the bluetooth address

Implements [AttysCommBase](#).

#### 4.1.3.5 start()

```
virtual void AttysComm::start ( ) [virtual]
```

Starts the data acquisition by starting the main thread. and sending possibly init commands.

Reimplemented from [AttysCommBase](#).


The documentation for this class was generated from the following file:

- AttysComm.h

## 4.2 AttysCommBase Class Reference

```
#include <AttysCommBase.h>
```

Inheritance diagram for AttysCommBase:





## Public Member Functions

- void [setAdc\\_samplingrate\\_index](#) (int idx)
- int [getSamplingRateInHz](#) ()
- int [getAdc\\_samplingrate\\_index](#) ()
- float [getADCFullScaleRange](#) (int channel)
- void [setAdc0\\_gain\\_index](#) (int idx)
- void [setAdc1\\_gain\\_index](#) (int idx)
- void [setBiasCurrent](#) (int currIndex)
- int [getBiasCurrent](#) ()
- void [enableCurrents](#) (int pos\_ch1, int neg\_ch1, int pos\_ch2)
- float [getAccelFullScaleRange](#) ()
- void [setAccel\\_full\\_scale\\_index](#) (int idx)
- float [getMagFullScaleRange](#) ()
- int [getIsCharging](#) ()
- virtual void [connect](#) ()=0
- virtual void [start](#) ()
- virtual void [closeSocket](#) ()=0
- int [hasActiveConnection](#) ()
- sample\_p [getSampleFromBuffer](#) ()
- int [hasSampleAvailable](#) ()
- void [resetRingbuffer](#) ()
- void [registerCallback](#) (AttysCommListener \*f)
- void [unregisterCallback](#) ()
- void [registerMessageCallback](#) (AttysCommMessage \*f)
- void [unregisterMessageCallback](#) ()
- void [quit](#) ()
- virtual unsigned char \* [getBluetoothBinaryAddress](#) ()=0
- virtual void [getBluetoothAddressString](#) (char \*s)=0

## Static Public Member Functions

- static float [phys2temperature](#) (float volt)

## Public Attributes

- const std::string [CHANNEL\\_DESCRIPTION](#) [NCHANNELS]
- const std::string [CHANNEL\\_SHORT\\_DESCRIPTION](#) [NCHANNELS]
- std::string const [CHANNEL\\_UNITS](#) [NCHANNELS]
- const int [ADC\\_SAMPLINGRATE](#) [4] = { 125, 250, 500, 1000 }
- const int [ADC\\_GAIN\\_FACTOR](#) [7] = { 6, 1, 2, 3, 4, 8, 12 }
- const float [ADC\\_REF](#) = 2.42F
- const float [oneG](#) = 9.80665F
- const float [ACCEL\\_FULL\\_SCALE](#) [4] = { 2 \* [oneG](#), 4 \* [oneG](#), 8 \* [oneG](#), 16 \* [oneG](#) }
- const float [MAG\\_FULL\\_SCALE](#) = 4800.0E-6F

## Static Public Attributes

- static const int [NCHANNELS](#) = 10
- static const int [nMem](#) = 1000 \* 10
- static const int [INDEX\\_Acceleration\\_X](#) = 0
- static const int [INDEX\\_Acceleration\\_Y](#) = 1
- static const int [INDEX\\_Acceleration\\_Z](#) = 2
- static const int [INDEX\\_Magnetic\\_field\\_X](#) = 3
- static const int [INDEX\\_Magnetic\\_field\\_Y](#) = 4
- static const int [INDEX\\_Magnetic\\_field\\_Z](#) = 5
- static const int [INDEX\\_Analogue\\_channel\\_1](#) = 6
- static const int [INDEX\\_Analogue\\_channel\\_2](#) = 7
- static const int [INDEX\\_GPIO0](#) = 8
- static const int [INDEX\\_GPIO1](#) = 9
- static const int [ADC\\_RATE\\_125HZ](#) = 0
- static const int [ADC\\_RATE\\_250HZ](#) = 1
- static const int [ADC\\_RATE\\_500Hz](#) = 2
- static const int [ADC\\_DEFAULT\\_RATE](#) = [ADC\\_RATE\\_250HZ](#)
- static const int [ADC\\_GAIN\\_6](#) = 0
- static const int [ADC\\_GAIN\\_1](#) = 1
- static const int [ADC\\_GAIN\\_2](#) = 2
- static const int [ADC\\_GAIN\\_3](#) = 3
- static const int [ADC\\_GAIN\\_4](#) = 4
- static const int [ADC\\_GAIN\\_8](#) = 5
- static const int [ADC\\_GAIN\\_12](#) = 6
- static const int [ADC\\_CURRENT\\_6NA](#) = 0
- static const int [ADC\\_CURRENT\\_22NA](#) = 1
- static const int [ADC\\_CURRENT\\_6UA](#) = 2
- static const int [ADC\\_CURRENT\\_22UA](#) = 3
- static const int [ADC\\_MUX\\_NORMAL](#) = 0
- static const int [ADC\\_MUX\\_SHORT](#) = 1
- static const int [ADC\\_MUX\\_SUPPLY](#) = 3
- static const int [ADC\\_MUX\\_TEMPERATURE](#) = 4
- static const int [ADC\\_MUX\\_TEST\\_SIGNAL](#) = 5
- static const int [ADC\\_MUX\\_ECG\\_EINTHOVEN](#) = 6
- static const int [ACCEL\\_2G](#) = 0
- static const int [ACCEL\\_4G](#) = 1
- static const int [ACCEL\\_8G](#) = 2
- static const int [ACCEL\\_16G](#) = 3
- static const int [MESSAGE\\_CONNECTED](#) = 0
- static const int [MESSAGE\\_ERROR](#) = 1
- static const int [MESSAGE\\_TIMEOUT](#) = 7
- static const int [MESSAGE\\_RECONNECTED](#) = 8
- static const int [MESSAGE\\_RECEIVING\\_DATA](#) = 9

## 4.2.1 Detailed Description

Platform independent definitions for the Attys

## 4.2.2 Member Function Documentation

#### 4.2.2.1 closeSocket()

```
virtual void AttysCommBase::closeSocket ( ) [pure virtual]
```

closes socket safely

Implemented in [AttysComm](#).

#### 4.2.2.2 connect()

```
virtual void AttysCommBase::connect ( ) [pure virtual]
```

connects to the Attys by opening the socket throws exception if it fails.

Implemented in [AttysComm](#).

#### 4.2.2.3 enableCurrents()

```
void AttysCommBase::enableCurrents (
    int pos_ch1,
    int neg_ch1,
    int pos_ch2 ) [inline]
```

Switches bias currents on

#### 4.2.2.4 getAccelFullScaleRange()

```
float AttysCommBase::getAccelFullScaleRange ( ) [inline]
```

Returns the accelerometer current full scale reading in m/s<sup>2</sup>

#### 4.2.2.5 getAdc\_samplingrate\_index()

```
int AttysCommBase::getAdc_samplingrate_index ( ) [inline]
```

Gets the sampling rate in form for the index.

#### 4.2.2.6 getADCFullScaleRange()

```
float AttysCommBase::getADCFullScaleRange (
    int channel ) [inline]
```

Gets the ADC full range. This depends on the gain setting of the ADC.

#### 4.2.2.7 getBiasCurrent()

```
int AttysCommBase::getBiasCurrent ( ) [inline]
```

Gets the bias current as in index.

#### 4.2.2.8 getBluetoothAdressString()

```
virtual void AttysCommBase::getBluetoothAdressString (
    char * s ) [pure virtual]
```

returns the Mac address as a string

Implemented in [AttysComm](#).

#### 4.2.2.9 getBluetoothBinaryAdress()

```
virtual unsigned char* AttysCommBase::getBluetoothBinaryAdress ( ) [pure virtual]
```

returns an array of 14 bytes of the bluetooth address

Implemented in [AttysComm](#).

#### 4.2.2.10 getIsCharging()

```
int AttysCommBase::getIsCharging ( ) [inline]
```

Charging indicator. Returns one if charging.

#### 4.2.2.11 getMagFullScaleRange()

```
float AttysCommBase::getMagFullScaleRange ( ) [inline]
```

Returns the full scale magnetometer in Tesla

#### 4.2.2.12 getSampleFromBuffer()

```
sample_p AttysCommBase::getSampleFromBuffer ( )
```

Gets a sample from the ringbuffer. This is a C array of all samples.

#### 4.2.2.13 getSamplingRateInHz()

```
int AttysCommBase::getSamplingRateInHz ( ) [inline]
```

Gets the sampling rate in Hz (not index number)

#### 4.2.2.14 hasActiveConnection()

```
int AttysCommBase::hasActiveConnection ( ) [inline]
```

Returns 1 if the connection is active.

#### 4.2.2.15 hasSampleAvailable()

```
int AttysCommBase::hasSampleAvailable ( ) [inline]
```

Is one if samples are available in the ringbuffer

#### 4.2.2.16 phys2temperature()

```
static float AttysCommBase::phys2temperature (
    float volt ) [inline], [static]
```

Temperature

#### 4.2.2.17 quit()

```
void AttysCommBase::quit ( )
```

Call this from the main activity to shutdown the connection

#### 4.2.2.18 registerCallback()

```
void AttysCommBase::registerCallback (
    AttysCommListener * f ) [inline]
```

Realtime callback function which is called whenever a sample has arrived. Implemented as an interface.

#### 4.2.2.19 registerMessageCallback()

```
void AttysCommBase::registerMessageCallback (
    AttysCommMessage * f ) [inline]
```

Callback function which is called whenever a special error/event has occurred.

#### 4.2.2.20 resetRingbuffer()

```
void AttysCommBase::resetRingbuffer ( ) [inline]
```

Resets the ringbuffer to zero content

#### 4.2.2.21 setAccel\_full\_scale\_index()

```
void AttysCommBase::setAccel_full_scale_index (
    int idx ) [inline]
```

Sets the accelerometer full scale range using the index.

#### 4.2.2.22 setAdc0\_gain\_index()

```
void AttysCommBase::setAdc0_gain_index (
    int idx ) [inline]
```

Gets the gain index for ADC1

**4.2.2.23 setAdc1\_gain\_index()**

```
void AttysCommBase::setAdc1_gain_index (
    int idx ) [inline]
```

Gets the gain index for ADC2

**4.2.2.24 setAdc\_samplingrate\_index()**

```
void AttysCommBase::setAdc_samplingrate_index (
    int idx ) [inline]
```

Sets the sampling rate using the sampling rate index numbers.

**4.2.2.25 setBiasCurrent()**

```
void AttysCommBase::setBiasCurrent (
    int currIndex ) [inline]
```

Sets the bias current which can be switched on.

**4.2.2.26 start()**

```
virtual void AttysCommBase::start ( ) [inline], [virtual]
```

Starts the data acquisition by starting the main thread. and sending possibly init commands.

Reimplemented in [AttysComm](#).

**4.2.2.27 unregisterCallback()**

```
void AttysCommBase::unregisterCallback ( ) [inline]
```

Unregister the sample callback

**4.2.2.28 unregisterMessageCallback()**

```
void AttysCommBase::unregisterMessageCallback ( ) [inline]
```

Unregister the message callback.

**4.2.3 Member Data Documentation****4.2.3.1 ACCEL\_16G**

```
const int AttysCommBase::ACCEL_16G = 3 [static]
```

Setting full scale range of the accelerometer to 16G

#### 4.2.3.2 ACCEL\_2G

```
const int AttysCommBase::ACCEL_2G = 0 [static]
```

Setting full scale range of the accelerometer to 2G

#### 4.2.3.3 ACCEL\_4G

```
const int AttysCommBase::ACCEL_4G = 1 [static]
```

Setting full scale range of the accelerometer to 4G

#### 4.2.3.4 ACCEL\_8G

```
const int AttysCommBase::ACCEL_8G = 2 [static]
```

Setting full scale range of the accelerometer to 8G

#### 4.2.3.5 ACCEL\_FULL\_SCALE

```
const float AttysCommBase::ACCEL_FULL_SCALE[4] = { 2 * oneG, 4 * oneG, 8 * oneG, 16 * oneG }
```

Mapping of the index to the full scale accelerations

#### 4.2.3.6 ADC\_CURRENT\_22NA

```
const int AttysCommBase::ADC_CURRENT_22NA = 1 [static]
```

Bias current of 22nA

#### 4.2.3.7 ADC\_CURRENT\_22UA

```
const int AttysCommBase::ADC_CURRENT_22UA = 3 [static]
```

Bias current of 22uA

#### 4.2.3.8 ADC\_CURRENT\_6NA

```
const int AttysCommBase::ADC_CURRENT_6NA = 0 [static]
```

Bias current of 6nA

#### 4.2.3.9 ADC\_CURRENT\_6UA

```
const int AttysCommBase::ADC_CURRENT_6UA = 2 [static]
```

Bias current of 6uA

#### 4.2.3.10 ADC\_DEFAULT\_RATE

```
const int AttysCommBase::ADC_DEFAULT_RATE = ADC_RATE_250HZ [static]
```

Constant defining the default sampling rate (250Hz)

#### 4.2.3.11 ADC\_GAIN\_1

```
const int AttysCommBase::ADC_GAIN_1 = 1 [static]
```

Gain index setting it to gain 6.

#### 4.2.3.12 ADC\_GAIN\_12

```
const int AttysCommBase::ADC_GAIN_12 = 6 [static]
```

Gain index setting it to gain 6.

#### 4.2.3.13 ADC\_GAIN\_2

```
const int AttysCommBase::ADC_GAIN_2 = 2 [static]
```

Gain index setting it to gain 2.

#### 4.2.3.14 ADC\_GAIN\_3

```
const int AttysCommBase::ADC_GAIN_3 = 3 [static]
```

Gain index setting it to gain 3.

#### 4.2.3.15 ADC\_GAIN\_4

```
const int AttysCommBase::ADC_GAIN_4 = 4 [static]
```

Gain index setting it to gain 4.

#### 4.2.3.16 ADC\_GAIN\_6

```
const int AttysCommBase::ADC_GAIN_6 = 0 [static]
```

Gain index setting it to gain 6.

#### 4.2.3.17 ADC\_GAIN\_8

```
const int AttysCommBase::ADC_GAIN_8 = 5 [static]
```

Gain index setting it to gain 5.

#### 4.2.3.18 ADC\_GAIN\_FACTOR

```
const int AttysCommBase::ADC_GAIN_FACTOR[7] = { 6, 1, 2, 3, 4, 8, 12 }
```

Mmapping between index and actual gain.

#### 4.2.3.19 ADC\_MUX\_ECG\_EINTHOVEN

```
const int AttysCommBase::ADC_MUX_ECG_EINTHOVEN = 6 [static]
```

Multiplexer routing: both positive ADC inputs are connected together



#### 4.2.3.20 ADC\_MUX\_NORMAL

```
const int AttysCommBase::ADC_MUX_NORMAL = 0 [static]
```

Multiplexer routing is normal: ADC1 and ADC2 are connected to the sigma/delta

#### 4.2.3.21 ADC\_MUX\_SHORT

```
const int AttysCommBase::ADC_MUX_SHORT = 1 [static]
```

Multiplexer routing: inputs are short circuited

#### 4.2.3.22 ADC\_MUX\_SUPPLY

```
const int AttysCommBase::ADC_MUX_SUPPLY = 3 [static]
```

Multiplexer routing: inputs are connected to power supply

#### 4.2.3.23 ADC\_MUX\_TEMPERATURE

```
const int AttysCommBase::ADC_MUX_TEMPERATURE = 4 [static]
```

Multiplexer routing: ADC measures internal temperature

#### 4.2.3.24 ADC\_MUX\_TEST\_SIGNAL

```
const int AttysCommBase::ADC_MUX_TEST_SIGNAL = 5 [static]
```

Multiplexer routing: ADC measures test signal

#### 4.2.3.25 ADC\_RATE\_125HZ

```
const int AttysCommBase::ADC_RATE_125HZ = 0 [static]
```

Constant defining sampling rate of 125Hz

#### 4.2.3.26 ADC\_RATE\_250HZ

```
const int AttysCommBase::ADC_RATE_250HZ = 1 [static]
```

Constant defining sampling rate of 250Hz

#### 4.2.3.27 ADC\_RATE\_500Hz

```
const int AttysCommBase::ADC_RATE_500Hz = 2 [static]
```

Constant defining sampling rate of 500Hz

#### 4.2.3.28 ADC\_REF

```
const float AttysCommBase::ADC_REF = 2.42F
```

The voltage reference of the ADC in volts

## 4.2.3.29 ADC\_SAMPLINGRATE

```
const int AttysCommBase::ADC_SAMPLINGRATE[4] = { 125, 250, 500, 1000 }
```

Array of the sampling rates mapping the index to the actual sampling rate.

## 4.2.3.30 CHANNEL\_DESCRIPTION

```
const std::string AttysCommBase::CHANNEL_DESCRIPTION[NCHANNELS]
```

**Initial value:**

```
= {
    "Acceleration X",
    "Acceleration Y",
    "Acceleration Z",
    "Magnetic field X",
    "Magnetic field Y",
    "Magnetic field Z",
    "Analogue channel 1",
    "Analogue channel 2",
    "DIN channel 0",
    "DIN channel 1",
    "Charging status"
}
```

Long descriptions of the channels in text form

## 4.2.3.31 CHANNEL\_SHORT\_DESCRIPTION

```
const std::string AttysCommBase::CHANNEL_SHORT_DESCRIPTION[NCHANNELS]
```

**Initial value:**

```
= {
    "Acc X",
    "Acc Y",
    "Acc Z",
    "Mag X",
    "Mag Y",
    "Mag Z",
    "ADC 1",
    "ADC 2",
    "DIN 0",
    "DIN 1",
}
```

Short descriptions of the channels in text form

## 4.2.3.32 CHANNEL\_UNITS

```
std::string const AttysCommBase::CHANNEL_UNITS[NCHANNELS]
```

**Initial value:**

```
= {
    "m/s^2",
    "m/s^2",
    "m/s^2",
    "T",
    "T",
    "T",
    "T",
    "V",
    "V",
    " ",
    " "
}
```

Units of the channels

#### 4.2.3.33 INDEX\_Acceleration\_X

```
const int AttysCommBase::INDEX_Acceleration_X = 0 [static]
```

Channel index for X Acceleration

#### 4.2.3.34 INDEX\_Acceleration\_Y

```
const int AttysCommBase::INDEX_Acceleration_Y = 1 [static]
```

Channel index for Y Acceleration

#### 4.2.3.35 INDEX\_Acceleration\_Z

```
const int AttysCommBase::INDEX_Acceleration_Z = 2 [static]
```

Channel index for Z Acceleration

#### 4.2.3.36 INDEX\_Analogue\_channel\_1

```
const int AttysCommBase::INDEX_Analogue_channel_1 = 6 [static]
```

Index of analogue channel 1

#### 4.2.3.37 INDEX\_Analogue\_channel\_2

```
const int AttysCommBase::INDEX_Analogue_channel_2 = 7 [static]
```

Index of analogue channel 2

#### 4.2.3.38 INDEX\_GPIO0

```
const int AttysCommBase::INDEX_GPIO0 = 8 [static]
```

Index of the internal GPIO pin 1

#### 4.2.3.39 INDEX\_GPIO1

```
const int AttysCommBase::INDEX_GPIO1 = 9 [static]
```

Index of the internal GPIO pin 2

#### 4.2.3.40 INDEX\_Magnetic\_field\_X

```
const int AttysCommBase::INDEX_Magnetic_field_X = 3 [static]
```

Magnetic field in X direction

#### 4.2.3.41 INDEX\_Magnetic\_field\_Y

```
const int AttysCommBase::INDEX_Magnetic_field_Y = 4 [static]
```

Magnetic field in Y direction

#### 4.2.3.42 INDEX\_Magnetic\_field\_Z

```
const int AttysCommBase::INDEX_Magnetic_field_Z = 5 [static]
```

Magnetic field in Z direction

#### 4.2.3.43 MAG\_FULL\_SCALE

```
const float AttysCommBase::MAG_FULL_SCALE = 4800.0E-6F
```

Full scale range of the magnetometer in Tesla

#### 4.2.3.44 MESSAGE\_CONNECTED

```
const int AttysCommBase::MESSAGE_CONNECTED = 0 [static]
```

Message callback: Connected.

#### 4.2.3.45 MESSAGE\_ERROR

```
const int AttysCommBase::MESSAGE_ERROR = 1 [static]
```

Message callback: Generic error.

#### 4.2.3.46 MESSAGE\_RECEIVING\_DATA

```
const int AttysCommBase::MESSAGE_RECEIVING_DATA = 9 [static]
```

Message callback: Receiving data.

#### 4.2.3.47 MESSAGE\_RECONNECTED

```
const int AttysCommBase::MESSAGE_RECONNECTED = 8 [static]
```

Message callback: Managed to reconnect.

#### 4.2.3.48 MESSAGE\_TIMEOUT

```
const int AttysCommBase::MESSAGE_TIMEOUT = 7 [static]
```

Message callback: Reception timeout detected by the watchdog.

#### 4.2.3.49 NCHANNELS

```
const int AttysCommBase::NCHANNELS = 10 [static]
```

Total number of channels per samples.

#### 4.2.3.50 nMem

```
const int AttysCommBase::nMem = 1000 * 10 [static]
```

Number of entries in the ringbuffer. Buffer for 10secs at 1kHz.

#### 4.2.3.51 oneG

```
const float AttysCommBase::oneG = 9.80665F
```

One g in m/s<sup>2</sup>

The documentation for this class was generated from the following file:

- AttysCommBase.h

### 4.3 AttysCommListener Struct Reference

```
#include <AttysCommBase.h>
```

#### 4.3.1 Detailed Description

callback when a sample has arrived

The documentation for this struct was generated from the following file:

- AttysCommBase.h

### 4.4 AttysCommMessage Struct Reference

```
#include <AttysCommBase.h>
```

#### 4.4.1 Detailed Description

callback when an error has occurred

The documentation for this struct was generated from the following file:

- AttysCommBase.h

## Index

ACCEL\_16G  
    AttysCommBase, 13  
ACCEL\_2G  
    AttysCommBase, 13  
ACCEL\_4G  
    AttysCommBase, 14  
ACCEL\_8G  
    AttysCommBase, 14  
ACCEL\_FULL\_SCALE  
    AttysCommBase, 14  
ADC\_CURRENT\_22NA  
    AttysCommBase, 14  
ADC\_CURRENT\_22UA  
    AttysCommBase, 14  
ADC\_CURRENT\_6NA  
    AttysCommBase, 14  
ADC\_CURRENT\_6UA  
    AttysCommBase, 14  
ADC\_DEFAULT\_RATE  
    AttysCommBase, 14  
ADC\_GAIN\_1  
    AttysCommBase, 14  
ADC\_GAIN\_12  
    AttysCommBase, 15  
ADC\_GAIN\_2  
    AttysCommBase, 15  
ADC\_GAIN\_3  
    AttysCommBase, 15  
ADC\_GAIN\_4  
    AttysCommBase, 15  
ADC\_GAIN\_6  
    AttysCommBase, 15  
ADC\_GAIN\_8  
    AttysCommBase, 15  
ADC\_GAIN\_FACTOR  
    AttysCommBase, 15  
ADC\_MUX\_ECG\_EINTHOVEN  
    AttysCommBase, 15  
ADC\_MUX\_NORMAL  
    AttysCommBase, 15  
ADC\_MUX\_SHORT  
    AttysCommBase, 16  
ADC\_MUX\_SUPPLY  
    AttysCommBase, 16  
ADC\_MUX\_TEMPERATURE  
    AttysCommBase, 16  
ADC\_MUX\_TEST\_SIGNAL  
    AttysCommBase, 16  
ADC\_RATE\_125HZ  
    AttysCommBase, 16  
ADC\_RATE\_250HZ  
    AttysCommBase, 16  
ADC\_RATE\_500Hz  
    AttysCommBase, 16  
ADC\_REF  
    AttysCommBase, 16  
ADC\_SAMPLINGRATE  
    AttysCommBase, 16  
AttysComm, 5  
    AttysComm, 6  
    closeSocket, 6  
    connect, 6  
    getBluetoothAdressString, 6  
    getBluetoothBinaryAdress, 7  
    start, 7  
AttysCommBase, 7  
    ACCEL\_16G, 13  
    ACCEL\_2G, 13  
    ACCEL\_4G, 14  
    ACCEL\_8G, 14  
    ACCEL\_FULL\_SCALE, 14  
    ADC\_CURRENT\_22NA, 14  
    ADC\_CURRENT\_22UA, 14  
    ADC\_CURRENT\_6NA, 14  
    ADC\_CURRENT\_6UA, 14  
    ADC\_DEFAULT\_RATE, 14  
    ADC\_GAIN\_1, 14  
    ADC\_GAIN\_12, 15  
    ADC\_GAIN\_2, 15  
    ADC\_GAIN\_3, 15  
    ADC\_GAIN\_4, 15  
    ADC\_GAIN\_6, 15  
    ADC\_GAIN\_8, 15  
    ADC\_GAIN\_FACTOR, 15  
    ADC\_MUX\_ECG\_EINTHOVEN, 15  
    ADC\_MUX\_NORMAL, 15  
    ADC\_MUX\_SHORT, 16  
    ADC\_MUX\_SUPPLY, 16  
    ADC\_MUX\_TEMPERATURE, 16  
    ADC\_MUX\_TEST\_SIGNAL, 16  
    ADC\_RATE\_125HZ, 16  
    ADC\_RATE\_250HZ, 16  
    ADC\_RATE\_500Hz, 16  
    ADC\_REF, 16  
    ADC\_SAMPLINGRATE, 16  
    CHANNEL\_DESCRIPTION, 17  
    CHANNEL\_SHORT\_DESCRIPTION, 17  
    CHANNEL\_UNITS, 17  
    closeSocket, 9  
    connect, 10  
    enableCurrents, 10  
    getADCFullScaleRange, 10  
    getAccelFullScaleRange, 10  
    getAdc\_samplingrate\_index, 10  
    getBiasCurrent, 10  
    getBluetoothAdressString, 10  
    getBluetoothBinaryAdress, 11  
    getIsCharging, 11  
    getMagFullScaleRange, 11  
    getSampleFromBuffer, 11

- getSamplingRateInHz, 11
- hasActiveConnection, 11
- hasSampleAvailable, 11
- INDEX\_Acceleration\_X, 17
- INDEX\_Acceleration\_Y, 18
- INDEX\_Acceleration\_Z, 18
- INDEX\_Analogue\_channel\_1, 18
- INDEX\_Analogue\_channel\_2, 18
- INDEX\_GPIO0, 18
- INDEX\_GPIO1, 18
- INDEX\_Magnetic\_field\_X, 18
- INDEX\_Magnetic\_field\_Y, 18
- INDEX\_Magnetic\_field\_Z, 18
- MAG\_FULL\_SCALE, 19
- MESSAGE\_CONNECTED, 19
- MESSAGE\_ERROR, 19
- MESSAGE\_RECEIVING\_DATA, 19
- MESSAGE\_RECONNECTED, 19
- MESSAGE\_TIMEOUT, 19
- NCHANNELS, 19
- nMem, 19
- oneG, 19
- phys2temperature, 12
- quit, 12
- registerCallback, 12
- registerMessageCallback, 12
- resetRingbuffer, 12
- setAccel\_full\_scale\_index, 12
- setAdc0\_gain\_index, 12
- setAdc1\_gain\_index, 12
- setAdc\_samplingrate\_index, 13
- setBiasCurrent, 13
- start, 13
- unregisterCallback, 13
- unregisterMessageCallback, 13
- AttysCommListener, 20
- AttysCommMessage, 20
- CHANNEL\_DESCRIPTION
  - AttysCommBase, 17
- CHANNEL\_SHORT\_DESCRIPTION
  - AttysCommBase, 17
- CHANNEL\_UNITS
  - AttysCommBase, 17
- closeSocket
  - AttysComm, 6
  - AttysCommBase, 9
- connect
  - AttysComm, 6
  - AttysCommBase, 10
- enableCurrents
  - AttysCommBase, 10
- getADCFullScaleRange
  - AttysCommBase, 10
- getAccelFullScaleRange
  - AttysCommBase, 10
- getAdc\_samplingrate\_index
  - AttysCommBase, 10
- getBiasCurrent
  - AttysCommBase, 10
- getBluetoothAddressString
  - AttysComm, 6
  - AttysCommBase, 10
- getBluetoothBinaryAddress
  - AttysComm, 7
  - AttysCommBase, 11
- getIsCharging
  - AttysCommBase, 11
- getMagFullScaleRange
  - AttysCommBase, 11
- getSampleFromBuffer
  - AttysCommBase, 11
- getSamplingRateInHz
  - AttysCommBase, 11
- hasActiveConnection
  - AttysCommBase, 11
- hasSampleAvailable
  - AttysCommBase, 11
- INDEX\_Acceleration\_X
  - AttysCommBase, 17
- INDEX\_Acceleration\_Y
  - AttysCommBase, 18
- INDEX\_Acceleration\_Z
  - AttysCommBase, 18
- INDEX\_Analogue\_channel\_1
  - AttysCommBase, 18
- INDEX\_Analogue\_channel\_2
  - AttysCommBase, 18
- INDEX\_GPIO0
  - AttysCommBase, 18
- INDEX\_GPIO1
  - AttysCommBase, 18
- INDEX\_Magnetic\_field\_X
  - AttysCommBase, 18
- INDEX\_Magnetic\_field\_Y
  - AttysCommBase, 18
- INDEX\_Magnetic\_field\_Z
  - AttysCommBase, 18
- MAG\_FULL\_SCALE
  - AttysCommBase, 19
- MESSAGE\_CONNECTED
  - AttysCommBase, 19
- MESSAGE\_ERROR
  - AttysCommBase, 19
- MESSAGE\_RECEIVING\_DATA
  - AttysCommBase, 19
- MESSAGE\_RECONNECTED
  - AttysCommBase, 19
- MESSAGE\_TIMEOUT
  - AttysCommBase, 19
- NCHANNELS
  - AttysCommBase, 19

nMem  
    AttysCommBase, [19](#)

oneG  
    AttysCommBase, [19](#)

phys2temperature  
    AttysCommBase, [12](#)

quit  
    AttysCommBase, [12](#)

registerCallback  
    AttysCommBase, [12](#)

registerMessageCallback  
    AttysCommBase, [12](#)

resetRingbuffer  
    AttysCommBase, [12](#)

setAccel\_full\_scale\_index  
    AttysCommBase, [12](#)

setAdc0\_gain\_index  
    AttysCommBase, [12](#)

setAdc1\_gain\_index  
    AttysCommBase, [12](#)

setAdc\_samplingrate\_index  
    AttysCommBase, [13](#)

setBiasCurrent  
    AttysCommBase, [13](#)

start  
    AttysComm, [7](#)  
    AttysCommBase, [13](#)

unregisterCallback  
    AttysCommBase, [13](#)

unregisterMessageCallback  
    AttysCommBase, [13](#)