attys-comm

Generated by Doxygen 1.8.13

Contents

1	Hierarchical Index			2		
	1.1	Class I	Hierarchy	2		
2	Clas	Class Index				
	2.1	Class I	_ist	2		
3	Class Documentation					
	3.1	AttysCo	omm Class Reference	2		
		3.1.1	Detailed Description	3		
		3.1.2	Constructor & Destructor Documentation	3		
		3.1.3	Member Function Documentation	3		
	3.2	AttysCo	ommBase Class Reference	5		
		3.2.1	Detailed Description	7		
		3.2.2	Constructor & Destructor Documentation	7		
		3.2.3	Member Function Documentation	8		
		3.2.4	Member Data Documentation	12		
	3.3	AttysCo	ommListener Struct Reference	19		
		3.3.1	Detailed Description	19		
		3.3.2	Member Function Documentation	19		
	3.4	AttysCo	ommMessage Struct Reference	19		
		3.4.1	Detailed Description	20		
		3.4.2	Member Function Documentation	20		
3.5 AttysScan Class Reference		can Class Reference	20			
		3.5.1	Detailed Description	20		
		3.5.2	Member Function Documentation	21		
		3.5.3	Member Data Documentation	21		
	3.6	AttysSo	canListener Struct Reference	22		
		3.6.1	Detailed Description	22		
Inc	dex			23		

1 Hierarchical Index

1.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

AttysCommBase	5
AttysComm	2
AttysCommListener	19
AttysCommMessage	19
AttysScan	20
AttysScanListener	22

2 Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:


AttysComm	2
AttysCommBase	5
AttysCommListener	19
AttysCommMessage	19
AttysScan	20
AttysScanListener	22

3 Class Documentation

3.1 AttysComm Class Reference

#include <AttysComm.h>

Inheritance diagram for AttysComm:



Public Member Functions

- AttysComm (void *_btAddr=NULL, int _btAddrLen=0)
- virtual void connect ()
- virtual void closeSocket ()
- virtual void run ()
- virtual void quit ()
- virtual void sendSyncCommand (const char *message, int waitForOK)
- virtual void sendInit ()
- virtual void start ()
- virtual void receptionTimeout ()
- unsigned char * getBluetoothBinaryAdress ()
- void getBluetoothAdressString (char *s)

Additional Inherited Members

3.1.1 Detailed Description

AttysComm contains all the neccessary comms to talk to the Attys on Linux, Windows and Mac. AttysComm class contains the device specific definitions and implements the abstract classes of AttysCommBase. See AttysCommBase for the definitions there. Instances of this class are automatically created by AttysScan and the user can ignore definitions here. All relevant user functions are in AttysCommBase. Use this class only if you have a fixed bluetooth address (Linux/Win) or a fixed bluetooth device (Mac) and won't need to scan for a bluetooth device.

3.1.2 Constructor & Destructor Documentation

3.1.2.1 AttysComm()

Constructor: Win/Linux: takes the bluetooth device structure and its length as an argument. For Mac: just a pointer to the bluetooth device (needs typecast to *_btAddr) and provide no length.

3.1.3 Member Function Documentation

3.1.3.1 closeSocket()

```
virtual void AttysComm::closeSocket ( ) [virtual]
```

Closes the socket safely.

Implements AttysCommBase.

```
3.1.3.2 connect()
```

virtual void AttysComm::connect () [virtual]

Connects to the Attys by opening the socket. Throws char* exception if it fails.

Implements AttysCommBase.

```
3.1.3.3 getBluetoothAdressString()
```

```
void AttysComm::getBluetoothAdressString ( {\tt char} \ * \ s \ ) \quad {\tt [virtual]}
```

returns the MAC address as a string.

Implements AttysCommBase.

3.1.3.4 getBluetoothBinaryAdress()

```
unsigned char* AttysComm::getBluetoothBinaryAdress ( ) [virtual]
```

Returns an array of 14 bytes of the bluetooth address.

Implements AttysCommBase.

```
3.1.3.5 quit()
```

```
virtual void AttysComm::quit ( ) [inline], [virtual]
```

Call this from the main activity to shut down the connection.

3.1.3.6 receptionTimeout()

```
virtual void AttysComm::receptionTimeout ( ) [virtual]
```

Called from the watchdog after a timeout. Do not call this directly.

Implements AttysCommBase.

```
3.1.3.7 run()
```

```
virtual void AttysComm::run ( ) [virtual]
```

Thread which does the data acquisition. Do not call directly.

Implements AttysCommBase.

3.1.3.8 sendInit()

```
virtual void AttysComm::sendInit ( ) [virtual]
```

Sends the init sequence to the Attys. Do not use unless you know exactly what you are doing.

Implements AttysCommBase.

3.1.3.9 sendSyncCommand()

Sends a command to the Attys. Do not use unless you know exactly what you are doing.

Implements AttysCommBase.

3.1.3.10 start()

```
virtual void AttysComm::start ( ) [virtual]
```

Starts the data acquisition by starting the main thread. and sending possibly init commands.

Reimplemented from AttysCommBase.


The documentation for this class was generated from the following file:

· AttysComm.h

3.2 AttysCommBase Class Reference

```
#include <AttysCommBase.h>
```

Inheritance diagram for AttysCommBase:



Public Member Functions

- AttysCommBase ()
- virtual ∼AttysCommBase ()
- void setAdc_samplingrate_index (int idx)
- int getSamplingRateInHz ()
- int getAdc samplingrate index ()
- float getADCFullScaleRange (int channel)
- void setAdc0_gain_index (int idx)
- void setAdc1_gain_index (int idx)
- void setBiasCurrent (int currIndex)
- int getBiasCurrent ()
- void enableCurrents (int pos_ch1, int neg_ch1, int pos_ch2)
- float getAccelFullScaleRange ()
- void setAccel full scale index (int idx)
- · float getMagFullScaleRange ()
- int getIsCharging ()
- virtual void connect ()=0
- · virtual void start ()
- virtual void closeSocket ()=0
- int hasActiveConnection ()
- sample_p getSampleFromBuffer ()
- int hasSampleAvailable ()
- void resetRingbuffer ()
- void registerCallback (AttysCommListener *f)
- void unregisterCallback ()
- void registerMessageCallback (AttysCommMessage *f)
- void unregisterMessageCallback ()
- · void quit ()
- virtual unsigned char * getBluetoothBinaryAdress ()=0
- virtual void getBluetoothAdressString (char *s)=0
- char * getAttysName ()
- void setAttysName (char *s)

Static Public Member Functions

• static float phys2temperature (float volt)

Public Attributes

- const std::string CHANNEL_DESCRIPTION [NCHANNELS]
- const std::string CHANNEL_SHORT_DESCRIPTION [NCHANNELS]
- std::string const CHANNEL UNITS [NCHANNELS]
- const int ADC_SAMPLINGRATE [4] = { 125, 250, 500, 1000 }
- const int ADC_GAIN_FACTOR [7] = { 6, 1, 2, 3, 4, 8, 12 }
- const float ADC REF = 2.42F
- const float oneG = 9.80665F
- const float ACCEL FULL SCALE [4] = { 2 * oneG, 4 * oneG, 8 * oneG, 16 * oneG }
- const float MAG_FULL_SCALE = 4800.0E-6F

Static Public Attributes

```
    static const int NCHANNELS = 10

static const int nMem = 1000 * 10

    static const int INDEX_Acceleration_X = 0

• static const int INDEX Acceleration Y = 1

    static const int INDEX Acceleration Z = 2

    static const int INDEX_Magnetic_field_X = 3

    static const int INDEX Magnetic field Y = 4

    static const int INDEX_Magnetic_field_Z = 5

• static const int INDEX_Analogue_channel_1 = 6

    static const int INDEX Analogue channel 2 = 7

    static const int INDEX GPIO0 = 8

    static const int INDEX_GPIO1 = 9

static const int ADC_RATE_125HZ = 0

    static const int ADC_RATE_250HZ = 1

    static const int ADC RATE 500Hz = 2

• static const int ADC_RATE_1000Hz = 3

    static const int ADC DEFAULT RATE = ADC RATE 250HZ

    static const int ADC GAIN 6 = 0

    static const int ADC GAIN 1 = 1

    static const int ADC_GAIN_2 = 2

• static const int ADC_GAIN_3 = 3

    static const int ADC GAIN 4 = 4

• static const int ADC GAIN 8 = 5

    static const int ADC GAIN 12 = 6

• static const int ADC_CURRENT_6NA = 0

    static const int ADC CURRENT 22NA = 1

    static const int ADC_CURRENT_6UA = 2

• static const int ADC_CURRENT_22UA = 3

    static const int ADC_MUX_NORMAL = 0

    static const int ADC MUX SHORT = 1

    static const int ADC_MUX_SUPPLY = 3

    static const int ADC_MUX_TEMPERATURE = 4

    static const int ADC_MUX_TEST_SIGNAL = 5

• static const int ADC MUX ECG EINTHOVEN = 6

    static const int ACCEL 2G = 0

    static const int ACCEL 4G = 1

    static const int ACCEL 8G = 2

    static const int ACCEL 16G = 3

    static const int MESSAGE_CONNECTED = 0

• static const int MESSAGE_ERROR = 1

    static const int MESSAGE TIMEOUT = 7
```

3.2.1 Detailed Description

Platform independent definitions for the Attys

static const int MESSAGE_RECONNECTED = 8
 static const int MESSAGE_RECEIVING_DATA = 9

3.2.2 Constructor & Destructor Documentation

3.2.2.1 AttysCommBase()

```
{\tt AttysCommBase::AttysCommBase \ (\ )}
```

Constructor which is overloaded by AttysComm.

3.2.2.2 ~AttysCommBase()

```
virtual AttysCommBase::~AttysCommBase ( ) [virtual]
```

Destructor which releases memory and closes any open connection.

3.2.3 Member Function Documentation

3.2.3.1 closeSocket()

```
virtual void AttysCommBase::closeSocket ( ) [pure virtual]
```

Closes the socket safely.

Implemented in AttysComm.

3.2.3.2 connect()

```
virtual void AttysCommBase::connect ( ) [pure virtual]
```

Connects to the Attys by opening the socket. Throws char* exception if it fails.

Implemented in AttysComm.

3.2.3.3 enableCurrents()

```
void AttysCommBase::enableCurrents (
    int pos_ch1,
    int neg_ch1,
    int pos_ch2 ) [inline]
```

Switches bias currents on

3.2.3.4 getAccelFullScaleRange()

```
float AttysCommBase::getAccelFullScaleRange ( ) [inline]
```

Returns the accelerometer current full scale reading in m/s $^{\wedge}$ 2.

```
3.2.3.5 getAdc_samplingrate_index()
```

```
int AttysCommBase::getAdc_samplingrate_index ( ) [inline]
```

Gets the sampling rate in form for the index.

3.2.3.6 getADCFullScaleRange()

Gets the ADC full range. This depends on the gain setting of the ADC.

3.2.3.7 getAttysName()

```
char* AttysCommBase::getAttysName ( ) [inline]
```

Returns the name of the Attys

3.2.3.8 getBiasCurrent()

```
int AttysCommBase::getBiasCurrent ( ) [inline]
```

Gets the bias current as in index.

3.2.3.9 getBluetoothAdressString()

```
virtual void AttysCommBase::getBluetoothAdressString ( {\tt char} \ * \ s \ ) \quad [{\tt pure} \ {\tt virtual}]
```

returns the MAC address as a string.

Implemented in AttysComm.

3.2.3.10 getBluetoothBinaryAdress()

```
virtual unsigned char* AttysCommBase::getBluetoothBinaryAdress ( ) [pure virtual]
```

Returns an array of 14 bytes of the bluetooth address.

Implemented in AttysComm.

3.2.3.11 getIsCharging()

```
int AttysCommBase::getIsCharging ( ) [inline]
```

Charging indicator. Returns one if charging.

```
3.2.3.12 getMagFullScaleRange()
```

```
float AttysCommBase::getMagFullScaleRange ( ) [inline]
```

Returns the full scale magnetometer in Tesla.

```
3.2.3.13 getSampleFromBuffer()
```

```
sample_p AttysCommBase::getSampleFromBuffer ( )
```

Gets a sample from the ringbuffer. This is a float* array of all channels.

3.2.3.14 getSamplingRateInHz()

```
int AttysCommBase::getSamplingRateInHz ( ) [inline]
```

Gets the sampling rate in Hz (not index number).

3.2.3.15 hasActiveConnection()

```
int AttysCommBase::hasActiveConnection ( ) [inline]
```

Returns one if the connection is active.

3.2.3.16 hasSampleAvailable()

```
int AttysCommBase::hasSampleAvailable ( ) [inline]
```

Is set to one if samples are available in the ringbuffer.

3.2.3.17 phys2temperature()

Temperature

3.2.3.18 quit()

```
void AttysCommBase::quit ( ) \,
```

Call this from the main activity to shut down the connection.

3.2.3.19 registerCallback()

Register a realtime callback function which is called whenever a sample has arrived. AttysCommListener is an abstract class which needs to implement hasSample().

3.2.3.20 registerMessageCallback()

Callback which is called whenever a special error/event has occurred.

```
3.2.3.21 resetRingbuffer()
```

```
void AttysCommBase::resetRingbuffer ( ) [inline]
```

Resets the ringbuffer to zero content.

```
3.2.3.22 setAccel_full_scale_index()
```

```
void AttysCommBase::setAccel_full_scale_index ( int \ \textit{idx} \ ) \ \ [inline]
```

Sets the accelerometer full scale range using the index.

```
3.2.3.23 setAdc0_gain_index()
```

Gets the gain index for ADC1.

```
3.2.3.24 setAdc1_gain_index()
```

Gets the gain index for ADC2.

3.2.3.25 setAdc_samplingrate_index()

Sets the sampling rate using the sampling rate index numbers.

3.2.3.26 setAttysName()

```
void AttysCommBase::setAttysName ( {\tt char} \ * \ s \ ) \quad [{\tt inline}]
```

Sets the name of the Attys

3.2.3.27 setBiasCurrent()

Sets the bias current which can be switched on.

```
3.2.3.28 start()
```

```
virtual void AttysCommBase::start ( ) [inline], [virtual]
```

Starts the data acquisition by starting the main thread. and sending possibly init commands.

Reimplemented in AttysComm.

3.2.3.29 unregisterCallback()

```
void AttysCommBase::unregisterCallback ( ) [inline]
```

Unregister the realtime sample callback.

3.2.3.30 unregisterMessageCallback()

```
void AttysCommBase::unregisterMessageCallback ( ) [inline]
```

Unregister the error/event callback.

3.2.4 Member Data Documentation

3.2.4.1 ACCEL_16G

```
const int AttysCommBase::ACCEL_16G = 3 [static]
```

Setting full scale range of the accelerometer to 16G.

3.2.4.2 ACCEL_2G

```
const int AttysCommBase::ACCEL_2G = 0 [static]
```

Setting full scale range of the accelerometer to 2G.

3.2.4.3 ACCEL_4G

```
const int AttysCommBase::ACCEL_4G = 1 [static]
```

Setting full scale range of the accelerometer to 4G.

3.2.4.4 ACCEL_8G

```
const int AttysCommBase::ACCEL_8G = 2 [static]
```

Setting full scale range of the accelerometer to 8G.

```
3.2.4.5 ACCEL_FULL_SCALE
const float AttysCommBase::ACCEL_FULL_SCALE[4] = { 2 * oneG, 4 * oneG, 8 * oneG, 16 * oneG }
Mapping of the index to the full scale accelerations.
3.2.4.6 ADC_CURRENT_22NA
const int AttysCommBase::ADC_CURRENT_22NA = 1 [static]
Bias current of 22nA.
3.2.4.7 ADC_CURRENT_22UA
const int AttysCommBase::ADC_CURRENT_22UA = 3 [static]
Bias current of 22uA.
3.2.4.8 ADC_CURRENT_6NA
const int AttysCommBase::ADC_CURRENT_6NA = 0 [static]
Bias current of 6nA.
3.2.4.9 ADC_CURRENT_6UA
const int AttysCommBase::ADC_CURRENT_6UA = 2 [static]
Bias current of 6uA.
3.2.4.10 ADC_DEFAULT_RATE
const int AttysCommBase::ADC_DEFAULT_RATE = ADC_RATE_250HZ [static]
Constant defining the default sampling rate (250Hz).
3.2.4.11 ADC_GAIN_1
const int AttysCommBase::ADC_GAIN_1 = 1 [static]
Gain index setting it to gain 6.
3.2.4.12 ADC_GAIN_12
const int AttysCommBase::ADC_GAIN_12 = 6 [static]
Gain index setting it to gain 6.
3.2.4.13 ADC_GAIN_2
const int AttysCommBase::ADC_GAIN_2 = 2 [static]
```

Gain index setting it to gain 2.

```
3.2.4.14 ADC_GAIN_3
```

```
const int AttysCommBase::ADC_GAIN_3 = 3 [static]
```

Gain index setting it to gain 3.

```
3.2.4.15 ADC_GAIN_4
```

```
const int AttysCommBase::ADC_GAIN_4 = 4 [static]
```

Gain index setting it to gain 4.

```
3.2.4.16 ADC_GAIN_6
```

```
const int AttysCommBase::ADC_GAIN_6 = 0 [static]
```

Gain index setting it to gain 6.

```
3.2.4.17 ADC_GAIN_8
```

```
const int AttysCommBase::ADC_GAIN_8 = 5 [static]
```

Gain index setting it to gain 5.

3.2.4.18 ADC_GAIN_FACTOR

```
const int AttysCommBase::ADC_GAIN_FACTOR[7] = { 6, 1, 2, 3, 4, 8, 12 }
```

Mmapping between index and actual gain.

3.2.4.19 ADC_MUX_ECG_EINTHOVEN

```
const int AttysCommBase::ADC_MUX_ECG_EINTHOVEN = 6 [static]
```

Muliplexer routing: both positive ADC inputs are connected together.

3.2.4.20 ADC_MUX_NORMAL

```
const int AttysCommBase::ADC_MUX_NORMAL = 0 [static]
```

Muliplexer routing is normal: ADC1 and ADC2 are connected to the sigma/delta.

3.2.4.21 ADC_MUX_SHORT

```
const int AttysCommBase::ADC_MUX_SHORT = 1 [static]
```

Muliplexer routing: inputs are short circuited.

3.2.4.22 ADC_MUX_SUPPLY

```
const int AttysCommBase::ADC_MUX_SUPPLY = 3 [static]
```

Muliplexer routing: inputs are connected to power supply.

3.2.4.23 ADC_MUX_TEMPERATURE

```
const int AttysCommBase::ADC_MUX_TEMPERATURE = 4 [static]
```

Muliplexer routing: ADC measures internal temperature.

3.2.4.24 ADC_MUX_TEST_SIGNAL

```
const int AttysCommBase::ADC_MUX_TEST_SIGNAL = 5 [static]
```

Muliplexer routing: ADC measures test signal.

3.2.4.25 ADC_RATE_1000Hz

```
const int AttysCommBase::ADC_RATE_1000Hz = 3 [static]
```

Constant defining sampling rate of 1000Hz (experimental).

3.2.4.26 ADC_RATE_125HZ

```
const int AttysCommBase::ADC_RATE_125HZ = 0 [static]
```

Constant defining sampling rate of 125Hz.

3.2.4.27 ADC_RATE_250HZ

```
const int AttysCommBase::ADC_RATE_250HZ = 1 [static]
```

Constant defining sampling rate of 250Hz.

3.2.4.28 ADC_RATE_500Hz

```
const int AttysCommBase::ADC_RATE_500Hz = 2 [static]
```

Constant defining sampling rate of 500Hz (experimental).

3.2.4.29 ADC_REF

```
const float AttysCommBase::ADC_REF = 2.42F
```

The voltage reference of the ADC in volts.

3.2.4.30 ADC_SAMPLINGRATE

```
const int AttysCommBase::ADC_SAMPLINGRATE[4] = { 125, 250, 500, 1000 }
```

Array of the sampling rates mapping the index to the actual sampling rate.

3.2.4.31 CHANNEL_DESCRIPTION

```
const std::string AttysCommBase::CHANNEL_DESCRIPTION[NCHANNELS]
```

Initial value:

```
= {
    "Acceleration X",
    "Acceleration Y",
    "Acceleration Z",
    "Magnetic field X",
    "Magnetic field Y",
    "Magnetic field Z",
    "Analogue channel 1",
    "Analogue channel 2"
    "DIN channel 0",
    "DIN channel 1",
    "Charging status"
```

Long descriptions of the channels.

3.2.4.32 CHANNEL_SHORT_DESCRIPTION

```
const std::string AttysCommBase::CHANNEL_SHORT_DESCRIPTION[NCHANNELS]
```

Initial value:

```
= {
    "Acc X",
    "Acc Y",
    "Acc X",
    "Mag X",
    "Mag Y",
    "Mag Z",
    "ADC 1",
    "ADC 2",
    "DIN 0",
    "DIN 1",
```

Short descriptions of the channels.

3.2.4.33 CHANNEL_UNITS

```
\verb|std::string| const AttysCommBase::CHANNEL_UNITS[NCHANNELS]|
```

Initial value:

```
= {
    "m/s^2",
    "m/s^2",
    "m/s^2",
    "T",
    "T",
    "T",
    "V",
    "V",
    "",
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """,
    """",
    """,
    """,
    """,
    """",
    """",
    """",
    """",
    """,
    """,
```

Units of the channels.

```
3.2.4.34 INDEX_Acceleration_X
const int AttysCommBase::INDEX_Acceleration_X = 0 [static]
Channel index for X Acceleration.
3.2.4.35 INDEX_Acceleration_Y
const int AttysCommBase::INDEX_Acceleration_Y = 1 [static]
Channel index for Y Acceleration.
3.2.4.36 INDEX_Acceleration_Z
const int AttysCommBase::INDEX_Acceleration_Z = 2 [static]
Channel index for Z Acceleration.
3.2.4.37 INDEX_Analogue_channel_1
const int AttysCommBase::INDEX_Analogue_channel_1 = 6 [static]
Index of analogue channel 1.
3.2.4.38 INDEX_Analogue_channel_2
const int AttysCommBase::INDEX_Analogue_channel_2 = 7 [static]
Index of analogue channel 2.
3.2.4.39 INDEX_GPIO0
const int AttysCommBase::INDEX_GPIO0 = 8 [static]
Index of the internal GPIO pin 1.
3.2.4.40 INDEX_GPIO1
const int AttysCommBase::INDEX_GPIO1 = 9 [static]
Index of the internal GPIO pin 2.
3.2.4.41 INDEX_Magnetic_field_X
const int AttysCommBase::INDEX_Magnetic_field_X = 3 [static]
Magnetic field in X direction.
3.2.4.42 INDEX_Magnetic_field_Y
const int AttysCommBase::INDEX_Magnetic_field_Y = 4 [static]
```

Magnetic field in Y direction.

```
3.2.4.43 INDEX_Magnetic_field_Z
```

```
const int AttysCommBase::INDEX_Magnetic_field_Z = 5 [static]
```

Magnetic field in Z direction.

3.2.4.44 MAG_FULL_SCALE

```
const float AttysCommBase::MAG_FULL_SCALE = 4800.0E-6F
```

Full scale range of the magnetometer in Tesla.

3.2.4.45 MESSAGE CONNECTED

```
const int AttysCommBase::MESSAGE_CONNECTED = 0 [static]
```

Message callback: Connected.

3.2.4.46 MESSAGE_ERROR

```
const int AttysCommBase::MESSAGE_ERROR = 1 [static]
```

Message callback: Generic error.

3.2.4.47 MESSAGE_RECEIVING_DATA

```
const int AttysCommBase::MESSAGE_RECEIVING_DATA = 9 [static]
```

Message callback: Receiving data. One off once the connection has been set up.

3.2.4.48 MESSAGE_RECONNECTED

```
const int AttysCommBase::MESSAGE_RECONNECTED = 8 [static]
```

Message callback: Managed to reconnect.

3.2.4.49 MESSAGE_TIMEOUT

```
const int AttysCommBase::MESSAGE_TIMEOUT = 7 [static]
```

Message callback: Reception timeout detected by the watchdog.

3.2.4.50 NCHANNELS

```
const int AttysCommBase::NCHANNELS = 10 [static]
```

Total number of channels per samples.

3.2.4.51 nMem

```
const int AttysCommBase::nMem = 1000 * 10 [static]
```

Number of entries in the ringbuffer. Buffer for 10secs at 1kHz.

3.2.4.52 oneG

```
const float AttysCommBase::oneG = 9.80665F
```

One g in m/s 2 .

The documentation for this class was generated from the following file:

· AttysCommBase.h

3.3 AttysCommListener Struct Reference

```
#include <AttysCommBase.h>
```

Public Member Functions

• virtual void hasSample (double, sample_p)=0

3.3.1 Detailed Description

Callback after a sample has arrived. The main class can for example inherit class and implement hasSample.

3.3.2 Member Function Documentation

3.3.2.1 hasSample()

Provides the timestamp and an array of all channels. This is an abstract method and needs to be overloaded with a real method doing the work.

The documentation for this struct was generated from the following file:

· AttysCommBase.h

3.4 AttysCommMessage Struct Reference

```
#include <AttysCommBase.h>
```

Public Member Functions

virtual void hasMessage (int, const char *)=0

3.4.1 Detailed Description

Callback after an error has occurred. This callback is in particular useful after a broken connection has been re-established.

3.4.2 Member Function Documentation

3.4.2.1 hasMessage()

Provides the error number and a text message about the error.

The documentation for this struct was generated from the following file:

· AttysCommBase.h

3.5 AttysScan Class Reference

```
#include <AttysScan.h>
```

Public Member Functions

- int scan (int maxAttys=1)
- void registerCallback (AttysScanListener *f)
- void unregisterCallback ()
- AttysComm * getAttysComm (int i)
- char * getAttysName (int i)
- int getNAttysDevices ()

Static Public Attributes

- static const int SCAN CONNECTED = 0
- static const int SCAN SEARCHING = 1
- static const int SCAN_NODEV = 2
- static const int SCAN SOCKETERR = 3
- static const int SCAN_CONNECTING = 4
- static const int SCAN_CONNECTERR = 5
- static const int MAX_ATTYS_DEVS = 4

3.5.1 Detailed Description

Scans for Attys and creates instances of AttysComm for every detected/paired Attys. There is no need to create instances of AttysComm yourself. This is done by this class automatically.

3.5.2 Member Function Documentation

3.5.2.1 getAttysComm()

Obtains the pointer to a valid AttysComm class which has been successfully detected while scanning.

3.5.2.2 getAttysName()

Gets the Attys name as reported by the bluetooth manager

3.5.2.3 getNAttysDevices()

```
int AttysScan::getNAttysDevices ( ) [inline]
```

Returns the number of Attys devices

3.5.2.4 registerCallback()

```
\begin{tabular}{ll} \beg
```

Register callback which reports the scanning progress for example for a splash screen.

3.5.2.5 scan()

```
int AttysScan::scan (
    int maxAttys = 1 )
```

Scans for the specified number of devices and connects to them. By default only for one Attys. returns 0 on success

3.5.2.6 unregisterCallback()

```
void AttysScan::unregisterCallback ( ) [inline]
```

Unregisters the callback

3.5.3 Member Data Documentation

3.5.3.1 MAX_ATTYS_DEVS

```
const int AttysScan::MAX_ATTYS_DEVS = 4 [static]
```

Max number of Attys Devices

3.5.3.2 SCAN_CONNECTED

```
const int AttysScan::SCAN_CONNECTED = 0 [static]
```

Message index that the connection to an attys has been successful.

3.5.3.3 SCAN_CONNECTERR

```
const int AttysScan::SCAN_CONNECTERR = 5 [static]
```

Connection error during scanning

3.5.3.4 SCAN_CONNECTING

```
const int AttysScan::SCAN_CONNECTING = 4 [static]
```

In the process of connecting

3.5.3.5 SCAN_NODEV

```
const int AttysScan::SCAN_NODEV = 2 [static]
```

Message index that no Attys has been detected

3.5.3.6 SCAN_SEARCHING

```
const int AttysScan::SCAN_SEARCHING = 1 [static]
```

Message index that AttysScan is actively scanning

3.5.3.7 SCAN_SOCKETERR

```
const int AttysScan::SCAN_SOCKETERR = 3 [static]
```

Message that the socket could not be opened

The documentation for this class was generated from the following file:

· AttysScan.h

3.6 AttysScanListener Struct Reference

```
#include <AttysScan.h>
```

3.6.1 Detailed Description

Callback which reports the status of the scanner

The documentation for this struct was generated from the following file:

· AttysScan.h

Index

~AttysCommBase	ADC_RATE_250HZ
AttysCommBase, 8	AttysCommBase, 15
Atty500mmBase, 0	ADC RATE 500Hz
ACCEL_16G	AttysCommBase, 15
AttysCommBase, 12	ADC_REF
ACCEL 2G	AttysCommBase, 15
AttysCommBase, 12	ADC SAMPLINGRATE
ACCEL 4G	AttysCommBase, 15
AttysCommBase, 12	AttysComm, 2
ACCEL 8G	AttysComm, 3
AttysCommBase, 12	closeSocket, 3
ACCEL_FULL_SCALE	connect, 3
AttysCommBase, 12	getBluetoothAdressString, 4
ADC_CURRENT_22NA	getBluetoothBinaryAdress, 4
AttysCommBase, 13	quit, 4
ADC_CURRENT_22UA	receptionTimeout, 4
AttysCommBase, 13	run, 4
ADC CURRENT 6NA	sendInit, 4
AttysCommBase, 13	sendSyncCommand, 5
ADC_CURRENT_6UA	start, 5
AttysCommBase, 13	AttysCommBase, 5
ADC DEFAULT RATE	~AttysCommBase, 8
	-
AttysCommBase, 13 ADC GAIN 1	ACCEL_16G, 12
	ACCEL_2G, 12
AttysCommBase, 13	ACCEL 9C 12
ADC_GAIN_12	ACCEL_8G, 12
AttysCommBase, 13	ACCEL_FULL_SCALE, 12
ADC_GAIN_2	ADC_CURRENT_22NA, 13
AttysCommBase, 13	ADC_CURRENT_22UA, 13
ADC_GAIN_3	ADC_CURRENT_6NA, 13
AttysCommBase, 13	ADC_CURRENT_6UA, 13
ADC_GAIN_4	ADC_DEFAULT_RATE, 13
AttysCommBase, 14	ADC_GAIN_1, 13
ADC_GAIN_6	ADC_GAIN_12, 13
AttysCommBase, 14	ADC_GAIN_2, 13
ADC_GAIN_8	ADC_GAIN_3, 13
AttysCommBase, 14	ADC_GAIN_4, 14
ADC_GAIN_FACTOR	ADC_GAIN_6, 14
AttysCommBase, 14	ADC_GAIN_8, 14
ADC_MUX_ECG_EINTHOVEN	ADC_GAIN_FACTOR, 14
AttysCommBase, 14	ADC_MUX_ECG_EINTHOVEN, 14
ADC_MUX_NORMAL	ADC_MUX_NORMAL, 14
AttysCommBase, 14	ADC_MUX_SHORT, 14
ADC_MUX_SHORT	ADC_MUX_SUPPLY, 14
AttysCommBase, 14	ADC_MUX_TEMPERATURE, 14
ADC_MUX_SUPPLY	ADC_MUX_TEST_SIGNAL, 15
AttysCommBase, 14	ADC_RATE_1000Hz, 15
ADC_MUX_TEMPERATURE	ADC_RATE_125HZ, 15
AttysCommBase, 14	ADC_RATE_250HZ, 15
ADC_MUX_TEST_SIGNAL	ADC_RATE_500Hz, 15
AttysCommBase, 15	ADC_REF, 15
ADC_RATE_1000Hz	ADC_SAMPLINGRATE, 15
AttysCommBase, 15	AttysCommBase, 7
ADC_RATE_125HZ	CHANNEL_DESCRIPTION, 15
AttysCommBase, 15	CHANNEL_SHORT_DESCRIPTION, 16

24 INDEX

CHANNEL_UNITS, 16	MAX_ATTYS_DEVS, 21
closeSocket, 8	registerCallback, 21
connect, 8	SCAN_CONNECTERR, 22
enableCurrents, 8	SCAN_CONNECTED, 22
getADCFullScaleRange, 9	SCAN_CONNECTING, 22
getAccelFullScaleRange, 8	SCAN NODEV, 22
getAdc_samplingrate_index, 8	SCAN SEARCHING, 22
getAttysName, 9	SCAN_SOCKETERR, 22
getBiasCurrent, 9	scan, 21
getBluetoothAdressString, 9	unregisterCallback, 21
getBluetoothBinaryAdress, 9	AttysScanListener, 22
The state of the s	Allysocaricisterier, 22
getIsCharging, 9	CHANNEL DESCRIPTION
getMagFullScaleRange, 9	AttysCommBase, 15
getSampleFromBuffer, 10	CHANNEL SHORT DESCRIPTION
getSamplingRateInHz, 10	AttysCommBase, 16
hasActiveConnection, 10	CHANNEL UNITS
hasSampleAvailable, 10	-
INDEX_Acceleration_X, 16	AttysCommBase, 16
INDEX_Acceleration_Y, 17	closeSocket
INDEX_Acceleration_Z, 17	AttysComm, 3
INDEX_Analogue_channel_1, 17	AttysCommBase, 8
INDEX Analogue channel 2, 17	connect
INDEX_GPIO0, 17	AttysComm, 3
INDEX GPIO1, 17	AttysCommBase, 8
INDEX_Magnetic_field_X, 17	
INDEX_Magnetic_field_Y, 17	enableCurrents
INDEX Magnetic field Z, 17	AttysCommBase, 8
MAG FULL SCALE, 18	
MESSAGE_CONNECTED, 18	getADCFullScaleRange
MESSAGE_ERROR, 18	AttysCommBase, 9
	getAccelFullScaleRange
MESSAGE_RECEIVING_DATA, 18	AttysCommBase, 8
MESSAGE_RECONNECTED, 18	getAdc_samplingrate_index
MESSAGE_TIMEOUT, 18	AttysCommBase, 8
NCHANNELS, 18	getAttysComm
nMem, 18	AttysScan, 21
oneG, 18	getAttysName
phys2temperature, 10	AttysCommBase, 9
quit, 10	AttysScan, 21
registerCallback, 10	getBiasCurrent
registerMessageCallback, 10	AttysCommBase, 9
resetRingbuffer, 11	getBluetoothAdressString
setAccel_full_scale_index, 11	AttysComm, 4
setAdc0_gain_index, 11	AttysCommBase, 9
setAdc1_gain_index, 11	getBluetoothBinaryAdress
setAdc_samplingrate_index, 11	AttysComm, 4
setAttysName, 11	-
setBiasCurrent, 11	AttysCommBase, 9
start, 11	getIsCharging
unregisterCallback, 12	AttysCommBase, 9
unregisterGaliback, 12 unregisterMessageCallback, 12	getMagFullScaleRange
	AttysCommBase, 9
AttysCommListener, 19	getNAttysDevices
hasSample, 19	AttysScan, 21
AttysCommMessage, 19	getSampleFromBuffer
hasMessage, 20	AttysCommBase, 10
AttysScan, 20	getSamplingRateInHz
getAttysComm, 21	AttysCommBase, 10
getAttysName, 21	
getNAttysDevices, 21	hasActiveConnection

INDEX 25

A++ O P 40	A++ O = 4
AttysCommBase, 10	AttysComm, 4
hasMessage	registerCallback
AttysCommMessage, 20	AttysCommBase, 10
hasSample	AttysScan, 21
AttysCommListener, 19	registerMessageCallback
hasSampleAvailable	AttysCommBase, 10
AttysCommBase, 10	resetRingbuffer
	AttysCommBase, 11
INDEX_Acceleration_X	run
AttysCommBase, 16	AttysComm, 4
INDEX_Acceleration_Y	
AttysCommBase, 17	SCAN_CONNECTERR
INDEX_Acceleration_Z	AttysScan, 22
AttysCommBase, 17	SCAN_CONNECTED
INDEX_Analogue_channel_1	AttysScan, 22
AttysCommBase, 17	SCAN CONNECTING
INDEX_Analogue_channel_2	AttysScan, 22
AttysCommBase, 17	SCAN NODEV
INDEX GPIO0	AttysScan, 22
AttysCommBase, 17	SCAN_SEARCHING
INDEX GPIO1	AttysScan, 22
AttysCommBase, 17	SCAN_SOCKETERR
·	
INDEX_Magnetic_field_X	AttysScan, 22
AttysCommBase, 17	scan
INDEX_Magnetic_field_Y	AttysScan, 21
AttysCommBase, 17	sendInit
INDEX_Magnetic_field_Z	AttysComm, 4
AttysCommBase, 17	sendSyncCommand
MAG FULL COM F	AttysComm, 5
MAG_FULL_SCALE	setAccel_full_scale_index
AttysCommBase, 18	AttysCommBase, 11
MAX_ATTYS_DEVS	setAdc0_gain_index
AttysScan, 21	AttysCommBase, 11
MESSAGE_CONNECTED	setAdc1_gain_index
AttysCommBase, 18	AttysCommBase, 11
MESSAGE_ERROR	setAdc_samplingrate_index
AttysCommBase, 18	AttysCommBase, 11
MESSAGE_RECEIVING_DATA	setAttysName
AttysCommBase, 18	AttysCommBase, 11
MESSAGE RECONNECTED	setBiasCurrent
AttysCommBase, 18	AttysCommBase, 11
MESSAGE_TIMEOUT	start
AttysCommBase, 18	AttysComm, 5
, my occumination, 10	-
NCHANNELS	AttysCommBase, 11
AttysCommBase, 18	unregisterCallback
nMem	_
AttysCommBase, 18	AttysCommBase, 12
7 kky o commission, 10	AttysScan, 21
oneG	unregisterMessageCallback
AttysCommBase, 18	AttysCommBase, 12
, my cooning acc, 10	
phys2temperature	
AttysCommBase, 10	
quit	
AttysComm, 4	
AttysCommBase, 10	
•	
receptionTimeout	
•	