

attys-comm

Generated by Doxygen 1.8.13

Contents

1	Hierarchical Index	2
1.1	Class Hierarchy	2
2	Class Index	2
2.1	Class List	2
3	Class Documentation	2
3.1	AttysComm Class Reference	2
3.1.1	Detailed Description	3
3.1.2	Constructor & Destructor Documentation	3
3.1.3	Member Function Documentation	3
3.2	AttysCommBase Class Reference	5
3.2.1	Detailed Description	7
3.2.2	Constructor & Destructor Documentation	7
3.2.3	Member Function Documentation	8
3.2.4	Member Data Documentation	12
3.3	AttysCommListener Struct Reference	18
3.3.1	Detailed Description	19
3.3.2	Member Function Documentation	19
3.4	AttysCommMessage Struct Reference	19
3.4.1	Detailed Description	19
3.4.2	Member Function Documentation	19
3.5	AttysScan Class Reference	20
3.5.1	Detailed Description	20
3.5.2	Member Function Documentation	20
3.5.3	Member Data Documentation	21
3.6	AttysScanListener Struct Reference	22
3.6.1	Detailed Description	22
	Index	23

1 Hierarchical Index

1.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

AttysCommBase	5
AttysComm	2
AttysCommListener	18
AttysCommMessage	19
AttysScan	20
AttysScanListener	22

2 Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:


AttysComm	2
AttysCommBase	5
AttysCommListener	18
AttysCommMessage	19
AttysScan	20
AttysScanListener	22

3 Class Documentation

3.1 AttysComm Class Reference

```
#include <AttysComm.h>
```

Inheritance diagram for AttysComm:



Public Member Functions

- [AttysComm](#) (void * _btAddr=NULL, int _btAddrLen=0)
- virtual void [connect](#) ()
- virtual void [closeSocket](#) ()
- virtual void [run](#) ()
- virtual void [quit](#) ()
- virtual void [sendSyncCommand](#) (const char *message, int waitForOK)
- virtual void [sendInit](#) ()
- virtual void [start](#) ()
- virtual void [receptionTimeout](#) ()
- unsigned char * [getBluetoothBinaryAddress](#) ()
- void [getBluetoothAddressString](#) (char *s)

Additional Inherited Members

3.1.1 Detailed Description

[AttysComm](#) contains all the necessary comms to talk to the Attys on Linux, Windows and Mac.

[AttysComm](#) class contains the device specific definitions and implements the abstract classes of [AttysCommBase](#). See [AttysCommBase](#) for the definitions there. Instances of this class are automatically created by [AttysScan](#) and the user can ignore definitions here. All relevant user functions are in [AttysCommBase](#). Use this class only if you have a fixed bluetooth address (Linux/Win) or a fixed bluetooth device (Mac) and won't need to scan for a bluetooth device.

3.1.2 Constructor & Destructor Documentation

3.1.2.1 AttysComm()

```
AttysComm::AttysComm (
    void * _btAddr = NULL,
    int _btAddrLen = 0 ) [inline]
```

Constructor: Win/Linux: takes the bluetooth device structure and its length as an argument. For Mac: just a pointer to the bluetooth device (needs typecast to *_btAddr) and provide no length.

3.1.3 Member Function Documentation

3.1.3.1 closeSocket()

```
virtual void AttysComm::closeSocket ( ) [virtual]
```

Closes the socket safely.

Implements [AttysCommBase](#).

3.1.3.2 connect()

```
virtual void AttysComm::connect ( ) [virtual]
```

Connects to the Attys by opening the socket. Throws char* exception if it fails.

Implements [AttysCommBase](#).

3.1.3.3 getBluetoothAdressString()

```
void AttysComm::getBluetoothAdressString (
    char * s ) [virtual]
```

returns the MAC address as a string.

Implements [AttysCommBase](#).

3.1.3.4 getBluetoothBinaryAdress()

```
unsigned char* AttysComm::getBluetoothBinaryAdress ( ) [virtual]
```

Returns an array of 14 bytes of the bluetooth address.

Implements [AttysCommBase](#).

3.1.3.5 quit()

```
virtual void AttysComm::quit ( ) [inline], [virtual]
```

Call this from the main activity to shut down the connection.

3.1.3.6 receptionTimeout()

```
virtual void AttysComm::receptionTimeout ( ) [virtual]
```

Called from the watchdog after a timeout. Do not call this directly.

Implements [AttysCommBase](#).

3.1.3.7 run()

```
virtual void AttysComm::run ( ) [virtual]
```

Thread which does the data acquisition. Do not call directly.

Implements [AttysCommBase](#).

3.1.3.8 sendInit()

```
virtual void AttysComm::sendInit ( ) [virtual]
```

Sends the init sequence to the Attys. Do not use unless you know exactly what you are doing.

Implements [AttysCommBase](#).

3.1.3.9 sendSyncCommand()

```
virtual void AttysComm::sendSyncCommand (
    const char * message,
    int waitForOK ) [virtual]
```

Sends a command to the Attys. Do not use unless you know exactly what you are doing.

Implements [AttysCommBase](#).

3.1.3.10 start()

```
virtual void AttysComm::start ( ) [virtual]
```

Starts the data acquisition by starting the main thread. and sending possibly init commands.

Reimplemented from [AttysCommBase](#).


The documentation for this class was generated from the following file:

- [AttysComm.h](#)

3.2 AttysCommBase Class Reference

```
#include <AttysCommBase.h>
```

Inheritance diagram for AttysCommBase:



Public Member Functions

- [AttysCommBase](#) ()
- virtual [~AttysCommBase](#) ()
- void [setAdc_samplingrate_index](#) (int idx)
- int [getSamplingRateInHz](#) ()
- int [getAdc_samplingrate_index](#) ()
- float [getADCFullScaleRange](#) (int channel)
- void [setAdc0_gain_index](#) (int idx)
- void [setAdc1_gain_index](#) (int idx)
- void [setBiasCurrent](#) (int currIndex)
- int [getBiasCurrent](#) ()
- void [enableCurrents](#) (int pos_ch1, int neg_ch1, int pos_ch2)
- float [getAccelFullScaleRange](#) ()
- void [setAccel_full_scale_index](#) (int idx)
- float [getMagFullScaleRange](#) ()
- int [getIsCharging](#) ()
- virtual void [connect](#) ()=0
- virtual void [start](#) ()
- virtual void [closeSocket](#) ()=0
- int [hasActiveConnection](#) ()
- sample_p [getSampleFromBuffer](#) ()
- int [hasSampleAvailable](#) ()
- void [resetRingbuffer](#) ()
- void [registerCallback](#) ([AttysCommListener](#) *f)
- void [unregisterCallback](#) ()
- void [registerMessageCallback](#) ([AttysCommMessage](#) *f)
- void [unregisterMessageCallback](#) ()
- void [quit](#) ()
- virtual unsigned char * [getBluetoothBinaryAddress](#) ()=0
- virtual void [getBluetoothAddressString](#) (char *s)=0

Static Public Member Functions

- static float [phys2temperature](#) (float volt)

Public Attributes

- const std::string [CHANNEL_DESCRIPTION](#) [[NCHANNELS](#)]
- const std::string [CHANNEL_SHORT_DESCRIPTION](#) [[NCHANNELS](#)]
- std::string const [CHANNEL_UNITS](#) [[NCHANNELS](#)]
- const int [ADC_SAMPLINGRATE](#) [4] = { 125, 250, 500, 1000 }
- const int [ADC_GAIN_FACTOR](#) [7] = { 6, 1, 2, 3, 4, 8, 12 }
- const float [ADC_REF](#) = 2.42F
- const float [oneG](#) = 9.80665F
- const float [ACCEL_FULL_SCALE](#) [4] = { 2 * [oneG](#), 4 * [oneG](#), 8 * [oneG](#), 16 * [oneG](#) }
- const float [MAG_FULL_SCALE](#) = 4800.0E-6F

Static Public Attributes

- static const int [NCHANNELS](#) = 10
- static const int [nMem](#) = 1000 * 10
- static const int [INDEX_Acceleration_X](#) = 0
- static const int [INDEX_Acceleration_Y](#) = 1
- static const int [INDEX_Acceleration_Z](#) = 2
- static const int [INDEX_Magnetic_field_X](#) = 3
- static const int [INDEX_Magnetic_field_Y](#) = 4
- static const int [INDEX_Magnetic_field_Z](#) = 5
- static const int [INDEX_Analogue_channel_1](#) = 6
- static const int [INDEX_Analogue_channel_2](#) = 7
- static const int [INDEX_GPIO0](#) = 8
- static const int [INDEX_GPIO1](#) = 9
- static const int [ADC_RATE_125HZ](#) = 0
- static const int [ADC_RATE_250HZ](#) = 1
- static const int [ADC_RATE_500Hz](#) = 2
- static const int [ADC_RATE_1000Hz](#) = 3
- static const int [ADC_DEFAULT_RATE](#) = [ADC_RATE_250HZ](#)
- static const int [ADC_GAIN_6](#) = 0
- static const int [ADC_GAIN_1](#) = 1
- static const int [ADC_GAIN_2](#) = 2
- static const int [ADC_GAIN_3](#) = 3
- static const int [ADC_GAIN_4](#) = 4
- static const int [ADC_GAIN_8](#) = 5
- static const int [ADC_GAIN_12](#) = 6
- static const int [ADC_CURRENT_6NA](#) = 0
- static const int [ADC_CURRENT_22NA](#) = 1
- static const int [ADC_CURRENT_6UA](#) = 2
- static const int [ADC_CURRENT_22UA](#) = 3
- static const int [ADC_MUX_NORMAL](#) = 0
- static const int [ADC_MUX_SHORT](#) = 1
- static const int [ADC_MUX_SUPPLY](#) = 3
- static const int [ADC_MUX_TEMPERATURE](#) = 4
- static const int [ADC_MUX_TEST_SIGNAL](#) = 5
- static const int [ADC_MUX_ECG_EINTHOVEN](#) = 6
- static const int [ACCEL_2G](#) = 0
- static const int [ACCEL_4G](#) = 1
- static const int [ACCEL_8G](#) = 2
- static const int [ACCEL_16G](#) = 3
- static const int [MESSAGE_CONNECTED](#) = 0
- static const int [MESSAGE_ERROR](#) = 1
- static const int [MESSAGE_TIMEOUT](#) = 7
- static const int [MESSAGE_RECONNECTED](#) = 8
- static const int [MESSAGE_RECEIVING_DATA](#) = 9

3.2.1 Detailed Description

Platform independent definitions for the Attys

3.2.2 Constructor & Destructor Documentation

3.2.2.1 AttysCommBase()

```
AttysCommBase::AttysCommBase ( )
```

Constructor which is overloaded by [AttysComm](#).

3.2.2.2 ~AttysCommBase()

```
virtual AttysCommBase::~~AttysCommBase ( ) [virtual]
```

Destructor which releases memory and closes any open connection.

3.2.3 Member Function Documentation

3.2.3.1 closeSocket()

```
virtual void AttysCommBase::closeSocket ( ) [pure virtual]
```

Closes the socket safely.

Implemented in [AttysComm](#).

3.2.3.2 connect()

```
virtual void AttysCommBase::connect ( ) [pure virtual]
```

Connects to the Attys by opening the socket. Throws char* exception if it fails.

Implemented in [AttysComm](#).

3.2.3.3 enableCurrents()

```
void AttysCommBase::enableCurrents (
    int pos_ch1,
    int neg_ch1,
    int pos_ch2 ) [inline]
```

Switches bias currents on

3.2.3.4 getAccelFullScaleRange()

```
float AttysCommBase::getAccelFullScaleRange ( ) [inline]
```

Returns the accelerometer current full scale reading in m/s².

3.2.3.5 getAdc_samplingrate_index()

```
int AttysCommBase::getAdc_samplingrate_index ( ) [inline]
```

Gets the sampling rate in form for the index.

3.2.3.6 getADCFullScaleRange()

```
float AttysCommBase::getADCFullScaleRange (
    int channel ) [inline]
```

Gets the ADC full range. This depends on the gain setting of the ADC.

3.2.3.7 getBiasCurrent()

```
int AttysCommBase::getBiasCurrent ( ) [inline]
```

Gets the bias current as in index.

3.2.3.8 getBluetoothAdressString()

```
virtual void AttysCommBase::getBluetoothAdressString (
    char * s ) [pure virtual]
```

returns the MAC address as a string.

Implemented in [AttysComm](#).

3.2.3.9 getBluetoothBinaryAdress()

```
virtual unsigned char* AttysCommBase::getBluetoothBinaryAdress ( ) [pure virtual]
```

Returns an array of 14 bytes of the bluetooth address.

Implemented in [AttysComm](#).

3.2.3.10 getIsCharging()

```
int AttysCommBase::getIsCharging ( ) [inline]
```

Charging indicator. Returns one if charging.

3.2.3.11 getMagFullScaleRange()

```
float AttysCommBase::getMagFullScaleRange ( ) [inline]
```

Returns the full scale magnetometer in Tesla.

3.2.3.12 `getSampleFromBuffer()`

```
sample_p AttysCommBase::getSampleFromBuffer ( )
```

Gets a sample from the ringbuffer. This is a float* array of all channels.

3.2.3.13 `getSamplingRateInHz()`

```
int AttysCommBase::getSamplingRateInHz ( ) [inline]
```

Gets the sampling rate in Hz (not index number).

3.2.3.14 `hasActiveConnection()`

```
int AttysCommBase::hasActiveConnection ( ) [inline]
```

Returns one if the connection is active.

3.2.3.15 `hasSampleAvailable()`

```
int AttysCommBase::hasSampleAvailable ( ) [inline]
```

Is set to one if samples are available in the ringbuffer.

3.2.3.16 `phys2temperature()`

```
static float AttysCommBase::phys2temperature (
    float volt ) [inline], [static]
```

Temperature

3.2.3.17 `quit()`

```
void AttysCommBase::quit ( )
```

Call this from the main activity to shut down the connection.

3.2.3.18 `registerCallback()`

```
void AttysCommBase::registerCallback (
    AttysCommListener * f ) [inline]
```

Register a realtime callback function which is called whenever a sample has arrived. [AttysCommListener](#) is an abstract class which needs to implement `hasSample()`.

3.2.3.19 `registerMessageCallback()`

```
void AttysCommBase::registerMessageCallback (
    AttysCommMessage * f ) [inline]
```

Callback which is called whenever a special error/event has occurred.

3.2.3.20 resetRingbuffer()

```
void AttysCommBase::resetRingbuffer ( ) [inline]
```

Resets the ringbuffer to zero content.

3.2.3.21 setAccel_full_scale_index()

```
void AttysCommBase::setAccel_full_scale_index (
    int idx ) [inline]
```

Sets the accelerometer full scale range using the index.

3.2.3.22 setAdc0_gain_index()

```
void AttysCommBase::setAdc0_gain_index (
    int idx ) [inline]
```

Gets the gain index for ADC1.

3.2.3.23 setAdc1_gain_index()

```
void AttysCommBase::setAdc1_gain_index (
    int idx ) [inline]
```

Gets the gain index for ADC2.

3.2.3.24 setAdc_samplingrate_index()

```
void AttysCommBase::setAdc_samplingrate_index (
    int idx ) [inline]
```

Sets the sampling rate using the sampling rate index numbers.

3.2.3.25 setBiasCurrent()

```
void AttysCommBase::setBiasCurrent (
    int currIndex ) [inline]
```

Sets the bias current which can be switched on.

3.2.3.26 start()

```
virtual void AttysCommBase::start ( ) [inline], [virtual]
```

Starts the data acquisition by starting the main thread. and sending possibly init commands.

Reimplemented in [AttysComm](#).

3.2.3.27 unregisterCallback()

```
void AttysCommBase::unregisterCallback ( ) [inline]
```

Unregister the realtime sample callback.

3.2.3.28 unregisterMessageCallback()

```
void AttysCommBase::unregisterMessageCallback ( ) [inline]
```

Unregister the error/event callback.

3.2.4 Member Data Documentation

3.2.4.1 ACCEL_16G

```
const int AttysCommBase::ACCEL_16G = 3 [static]
```

Setting full scale range of the accelerometer to 16G.

3.2.4.2 ACCEL_2G

```
const int AttysCommBase::ACCEL_2G = 0 [static]
```

Setting full scale range of the accelerometer to 2G.

3.2.4.3 ACCEL_4G

```
const int AttysCommBase::ACCEL_4G = 1 [static]
```

Setting full scale range of the accelerometer to 4G.

3.2.4.4 ACCEL_8G

```
const int AttysCommBase::ACCEL_8G = 2 [static]
```

Setting full scale range of the accelerometer to 8G.

3.2.4.5 ACCEL_FULL_SCALE

```
const float AttysCommBase::ACCEL_FULL_SCALE[4] = { 2 * oneG, 4 * oneG, 8 * oneG, 16 * oneG }
```

Mapping of the index to the full scale accelerations.

3.2.4.6 ADC_CURRENT_22NA

```
const int AttysCommBase::ADC_CURRENT_22NA = 1 [static]
```

Bias current of 22nA.

3.2.4.7 ADC_CURRENT_22UA

```
const int AttysCommBase::ADC_CURRENT_22UA = 3 [static]
```

Bias current of 22uA.

3.2.4.8 ADC_CURRENT_6NA

```
const int AttysCommBase::ADC_CURRENT_6NA = 0 [static]
```

Bias current of 6nA.

3.2.4.9 ADC_CURRENT_6UA

```
const int AttysCommBase::ADC_CURRENT_6UA = 2 [static]
```

Bias current of 6uA.

3.2.4.10 ADC_DEFAULT_RATE

```
const int AttysCommBase::ADC_DEFAULT_RATE = ADC_RATE_250HZ [static]
```

Constant defining the default sampling rate (250Hz).

3.2.4.11 ADC_GAIN_1

```
const int AttysCommBase::ADC_GAIN_1 = 1 [static]
```

Gain index setting it to gain 6.

3.2.4.12 ADC_GAIN_12

```
const int AttysCommBase::ADC_GAIN_12 = 6 [static]
```

Gain index setting it to gain 6.

3.2.4.13 ADC_GAIN_2

```
const int AttysCommBase::ADC_GAIN_2 = 2 [static]
```

Gain index setting it to gain 2.

3.2.4.14 ADC_GAIN_3

```
const int AttysCommBase::ADC_GAIN_3 = 3 [static]
```

Gain index setting it to gain 3.

3.2.4.15 ADC_GAIN_4

```
const int AttysCommBase::ADC_GAIN_4 = 4 [static]
```

Gain index setting it to gain 4.

3.2.4.16 ADC_GAIN_6

```
const int AttysCommBase::ADC_GAIN_6 = 0 [static]
```

Gain index setting it to gain 6.

3.2.4.17 ADC_GAIN_8

```
const int AttysCommBase::ADC_GAIN_8 = 5 [static]
```

Gain index setting it to gain 5.

3.2.4.18 ADC_GAIN_FACTOR

```
const int AttysCommBase::ADC_GAIN_FACTOR[7] = { 6, 1, 2, 3, 4, 8, 12 }
```

Mmapping between index and actual gain.

3.2.4.19 ADC_MUX_ECG_EINTHOVEN

```
const int AttysCommBase::ADC_MUX_ECG_EINTHOVEN = 6 [static]
```

Multiplexer routing: both positive ADC inputs are connected together.

3.2.4.20 ADC_MUX_NORMAL

```
const int AttysCommBase::ADC_MUX_NORMAL = 0 [static]
```

Multiplexer routing is normal: ADC1 and ADC2 are connected to the sigma/delta.

3.2.4.21 ADC_MUX_SHORT

```
const int AttysCommBase::ADC_MUX_SHORT = 1 [static]
```

Multiplexer routing: inputs are short circuited.

3.2.4.22 ADC_MUX_SUPPLY

```
const int AttysCommBase::ADC_MUX_SUPPLY = 3 [static]
```

Multiplexer routing: inputs are connected to power supply.

3.2.4.23 ADC_MUX_TEMPERATURE

```
const int AttysCommBase::ADC_MUX_TEMPERATURE = 4 [static]
```

Multiplexer routing: ADC measures internal temperature.

3.2.4.24 ADC_MUX_TEST_SIGNAL

```
const int AttysCommBase::ADC_MUX_TEST_SIGNAL = 5 [static]
```

Multiplexer routing: ADC measures test signal.

3.2.4.25 ADC_RATE_1000Hz

```
const int AttysCommBase::ADC_RATE_1000Hz = 3 [static]
```

Constant defining sampling rate of 1000Hz (experimental).

3.2.4.26 ADC_RATE_125HZ

```
const int AttysCommBase::ADC_RATE_125HZ = 0 [static]
```

Constant defining sampling rate of 125Hz.

3.2.4.27 ADC_RATE_250HZ

```
const int AttysCommBase::ADC_RATE_250HZ = 1 [static]
```

Constant defining sampling rate of 250Hz.

3.2.4.28 ADC_RATE_500Hz

```
const int AttysCommBase::ADC_RATE_500Hz = 2 [static]
```

Constant defining sampling rate of 500Hz (experimental).

3.2.4.29 ADC_REF

```
const float AttysCommBase::ADC_REF = 2.42F
```

The voltage reference of the ADC in volts.

3.2.4.30 ADC_SAMPLINGRATE

```
const int AttysCommBase::ADC_SAMPLINGRATE[4] = { 125, 250, 500, 1000 }
```

Array of the sampling rates mapping the index to the actual sampling rate.

3.2.4.31 CHANNEL_DESCRIPTION

```
const std::string AttysCommBase::CHANNEL_DESCRIPTION[NCHANNELS]
```

Initial value:

```
= {
    "Acceleration X",
    "Acceleration Y",
    "Acceleration Z",
    "Magnetic field X",
    "Magnetic field Y",
    "Magnetic field Z",
    "Analogue channel 1",
    "Analogue channel 2",
    "DIN channel 0",
    "DIN channel 1",
    "Charging status"
}
```

Long descriptions of the channels.

3.2.4.32 CHANNEL_SHORT_DESCRIPTION

```
const std::string AttysCommBase::CHANNEL_SHORT_DESCRIPTION[NCHANNELS]
```

Initial value:

```
= {
    "Acc X",
    "Acc Y",
    "Acc Z",
    "Mag X",
    "Mag Y",
    "Mag Z",
    "ADC 1",
    "ADC 2",
    "DIN 0",
    "DIN 1",
}
```

Short descriptions of the channels.

3.2.4.33 CHANNEL_UNITS

```
std::string const AttysCommBase::CHANNEL_UNITS[NCHANNELS]
```

Initial value:

```
= {
    "m/s^2",
    "m/s^2",
    "m/s^2",
    "T",
    "T",
    "T",
    "T",
    "V",
    "V",
    " ",
    " ",
}
```

Units of the channels.

3.2.4.34 INDEX_Acceleration_X

```
const int AttysCommBase::INDEX_Acceleration_X = 0 [static]
```

Channel index for X Acceleration.

3.2.4.35 INDEX_Acceleration_Y

```
const int AttysCommBase::INDEX_Acceleration_Y = 1 [static]
```

Channel index for Y Acceleration.

3.2.4.36 INDEX_Acceleration_Z

```
const int AttysCommBase::INDEX_Acceleration_Z = 2 [static]
```

Channel index for Z Acceleration.

3.2.4.37 INDEX_Analogue_channel_1

```
const int AttysCommBase::INDEX_Analogue_channel_1 = 6 [static]
```

Index of analogue channel 1.

3.2.4.38 INDEX_Analogue_channel_2

```
const int AttysCommBase::INDEX_Analogue_channel_2 = 7 [static]
```

Index of analogue channel 2.

3.2.4.39 INDEX_GPIO0

```
const int AttysCommBase::INDEX_GPIO0 = 8 [static]
```

Index of the internal GPIO pin 1.

3.2.4.40 INDEX_GPIO1

```
const int AttysCommBase::INDEX_GPIO1 = 9 [static]
```

Index of the internal GPIO pin 2.

3.2.4.41 INDEX_Magnetic_field_X

```
const int AttysCommBase::INDEX_Magnetic_field_X = 3 [static]
```

Magnetic field in X direction.

3.2.4.42 INDEX_Magnetic_field_Y

```
const int AttysCommBase::INDEX_Magnetic_field_Y = 4 [static]
```

Magnetic field in Y direction.

3.2.4.43 INDEX_Magnetic_field_Z

```
const int AttysCommBase::INDEX_Magnetic_field_Z = 5 [static]
```

Magnetic field in Z direction.

3.2.4.44 MAG_FULL_SCALE

```
const float AttysCommBase::MAG_FULL_SCALE = 4800.0E-6F
```

Full scale range of the magnetometer in Tesla.

3.2.4.45 MESSAGE_CONNECTED

```
const int AttysCommBase::MESSAGE_CONNECTED = 0 [static]
```

Message callback: Connected.

3.2.4.46 MESSAGE_ERROR

```
const int AttysCommBase::MESSAGE_ERROR = 1 [static]
```

Message callback: Generic error.

3.2.4.47 MESSAGE_RECEIVING_DATA

```
const int AttysCommBase::MESSAGE_RECEIVING_DATA = 9 [static]
```

Message callback: Receiving data.

3.2.4.48 MESSAGE_RECONNECTED

```
const int AttysCommBase::MESSAGE_RECONNECTED = 8 [static]
```

Message callback: Managed to reconnect.

3.2.4.49 MESSAGE_TIMEOUT

```
const int AttysCommBase::MESSAGE_TIMEOUT = 7 [static]
```

Message callback: Reception timeout detected by the watchdog.

3.2.4.50 NCHANNELS

```
const int AttysCommBase::NCHANNELS = 10 [static]
```

Total number of channels per samples.

3.2.4.51 nMem

```
const int AttysCommBase::nMem = 1000 * 10 [static]
```

Number of entries in the ringbuffer. Buffer for 10secs at 1kHz.

3.2.4.52 oneG

```
const float AttysCommBase::oneG = 9.80665F
```

One g in m/s².

The documentation for this class was generated from the following file:

- AttysCommBase.h

3.3 AttysCommListener Struct Reference

```
#include <AttysCommBase.h>
```

Public Member Functions

- virtual void [hasSample](#) (double, sample_p)=0

3.3.1 Detailed Description

Callback after a sample has arrived. The main class can for example inherit class and implement hasSample.

3.3.2 Member Function Documentation

3.3.2.1 hasSample()

```
virtual void AttysCommListener::hasSample (
    double ,
    sample_p ) [pure virtual]
```

Provides the timestamp and an array of all channels. This is an abstract method and needs to be overloaded with a real method doing the work.

The documentation for this struct was generated from the following file:

- AttysCommBase.h

3.4 AttysCommMessage Struct Reference

```
#include <AttysCommBase.h>
```

Public Member Functions

- virtual void [hasMessage](#) (int, const char *)=0

3.4.1 Detailed Description

Callback after an error has occurred. This callback is in particular useful after a broken connection has been re-established.

3.4.2 Member Function Documentation

3.4.2.1 hasMessage()

```
virtual void AttysCommMessage::hasMessage (
    int ,
    const char * ) [pure virtual]
```

Provides the error number and a text message about the error.

The documentation for this struct was generated from the following file:

- AttysCommBase.h

3.5 AttysScan Class Reference

```
#include <AttysScan.h>
```

Public Member Functions

- int [scan](#) (int maxAttys=1)
- void [registerCallback](#) (AttysScanListener *f)
- void [unregisterCallback](#) ()
- AttysComm * [getAttysComm](#) (int i)
- char * [getAttysName](#) (int i)
- int [getNAttysDevices](#) ()

Static Public Attributes

- static const int [SCAN_CONNECTED](#) = 0
- static const int [SCAN_SEARCHING](#) = 1
- static const int [SCAN_NODEV](#) = 2
- static const int [SCAN_SOCKETERR](#) = 3
- static const int [SCAN_CONNECTING](#) = 4
- static const int [SCAN_CONNECTERR](#) = 5
- static const int [MAX_ATTYS_DEVS](#) = 4

3.5.1 Detailed Description

Scans for Attys and creates instances of [AttysComm](#) for every detected/paired Attys. There is no need to create instances of [AttysComm](#) yourself. This is done by this class automatically.

3.5.2 Member Function Documentation

3.5.2.1 getAttysComm()

```
AttysComm* AttysScan::getAttysComm (
    int i ) [inline]
```

Obtains the pointer to a valid [AttysComm](#) class which has been successfully detected while scanning.

3.5.2.2 getAttysName()

```
char* AttysScan::getAttysName (
    int i ) [inline]
```

Gets the Attys name as reported by the bluetooth manager

3.5.2.3 getNAttysDevices()

```
int AttysScan::getNAttysDevices ( ) [inline]
```

Returns the number of Attys devices

3.5.2.4 registerCallback()

```
void AttysScan::registerCallback (
    AttysScanListener * f ) [inline]
```

Register callback which reports the scanning progress for example for a splash screen.

3.5.2.5 scan()

```
int AttysScan::scan (
    int maxAttys = 1 )
```

Scans for the specified number of devices and connects to them. By default only for one Attys. returns 0 on success

3.5.2.6 unregisterCallback()

```
void AttysScan::unregisterCallback ( ) [inline]
```

Unregisters the callback

3.5.3 Member Data Documentation

3.5.3.1 MAX_ATTYS_DEVS

```
const int AttysScan::MAX_ATTYS_DEVS = 4 [static]
```

Max number of Attys Devices

3.5.3.2 SCAN_CONNECTED

```
const int AttysScan::SCAN_CONNECTED = 0 [static]
```

Message index that the connection to an attys has been successful.

3.5.3.3 SCAN_CONNECTERR

```
const int AttysScan::SCAN_CONNECTERR = 5 [static]
```

Connection error during scanning

3.5.3.4 SCAN_CONNECTING

```
const int AttysScan::SCAN_CONNECTING = 4 [static]
```

In the process of connecting

3.5.3.5 SCAN_NODEV

```
const int AttysScan::SCAN_NODEV = 2 [static]
```

Message index that no Attys has been detected

3.5.3.6 SCAN_SEARCHING

```
const int AttysScan::SCAN_SEARCHING = 1 [static]
```

Message index that [AttysScan](#) is actively scanning

3.5.3.7 SCAN_SOCKETERR

```
const int AttysScan::SCAN_SOCKETERR = 3 [static]
```

Message that the socket could not be opened

The documentation for this class was generated from the following file:

- [AttysScan.h](#)

3.6 AttysScanListener Struct Reference

```
#include <AttysScan.h>
```

3.6.1 Detailed Description

Callback which reports the status of the scanner

The documentation for this struct was generated from the following file:

- [AttysScan.h](#)

Index

~AttysCommBase
 AttysCommBase, [8](#)

ACCEL_16G
 AttysCommBase, [12](#)

ACCEL_2G
 AttysCommBase, [12](#)

ACCEL_4G
 AttysCommBase, [12](#)

ACCEL_8G
 AttysCommBase, [12](#)

ACCEL_FULL_SCALE
 AttysCommBase, [12](#)

ADC_CURRENT_22NA
 AttysCommBase, [12](#)

ADC_CURRENT_22UA
 AttysCommBase, [12](#)

ADC_CURRENT_6NA
 AttysCommBase, [13](#)

ADC_CURRENT_6UA
 AttysCommBase, [13](#)

ADC_DEFAULT_RATE
 AttysCommBase, [13](#)

ADC_GAIN_1
 AttysCommBase, [13](#)

ADC_GAIN_12
 AttysCommBase, [13](#)

ADC_GAIN_2
 AttysCommBase, [13](#)

ADC_GAIN_3
 AttysCommBase, [13](#)

ADC_GAIN_4
 AttysCommBase, [13](#)

ADC_GAIN_6
 AttysCommBase, [13](#)

ADC_GAIN_8
 AttysCommBase, [14](#)

ADC_GAIN_FACTOR
 AttysCommBase, [14](#)

ADC_MUX_ECG_EINTHOVEN
 AttysCommBase, [14](#)

ADC_MUX_NORMAL
 AttysCommBase, [14](#)

ADC_MUX_SHORT
 AttysCommBase, [14](#)

ADC_MUX_SUPPLY
 AttysCommBase, [14](#)

ADC_MUX_TEMPERATURE
 AttysCommBase, [14](#)

ADC_MUX_TEST_SIGNAL
 AttysCommBase, [14](#)

ADC_RATE_1000Hz
 AttysCommBase, [14](#)

ADC_RATE_125HZ
 AttysCommBase, [15](#)

ADC_RATE_250HZ
 AttysCommBase, [15](#)

ADC_RATE_500Hz
 AttysCommBase, [15](#)

ADC_REF
 AttysCommBase, [15](#)

ADC_SAMPLINGRATE
 AttysCommBase, [15](#)

AttysComm, [2](#)
 AttysComm, [3](#)
 closeSocket, [3](#)
 connect, [3](#)
 getBluetoothAdressString, [4](#)
 getBluetoothBinaryAdress, [4](#)
 quit, [4](#)
 receptionTimeout, [4](#)
 run, [4](#)
 sendInit, [4](#)
 sendSyncCommand, [5](#)
 start, [5](#)

AttysCommBase, [5](#)
 ~AttysCommBase, [8](#)
 ACCEL_16G, [12](#)
 ACCEL_2G, [12](#)
 ACCEL_4G, [12](#)
 ACCEL_8G, [12](#)
 ACCEL_FULL_SCALE, [12](#)
 ADC_CURRENT_22NA, [12](#)
 ADC_CURRENT_22UA, [12](#)
 ADC_CURRENT_6NA, [13](#)
 ADC_CURRENT_6UA, [13](#)
 ADC_DEFAULT_RATE, [13](#)
 ADC_GAIN_1, [13](#)
 ADC_GAIN_12, [13](#)
 ADC_GAIN_2, [13](#)
 ADC_GAIN_3, [13](#)
 ADC_GAIN_4, [13](#)
 ADC_GAIN_6, [13](#)
 ADC_GAIN_8, [14](#)
 ADC_GAIN_FACTOR, [14](#)
 ADC_MUX_ECG_EINTHOVEN, [14](#)
 ADC_MUX_NORMAL, [14](#)
 ADC_MUX_SHORT, [14](#)
 ADC_MUX_SUPPLY, [14](#)
 ADC_MUX_TEMPERATURE, [14](#)
 ADC_MUX_TEST_SIGNAL, [14](#)
 ADC_RATE_1000Hz, [14](#)
 ADC_RATE_125HZ, [15](#)
 ADC_RATE_250HZ, [15](#)
 ADC_RATE_500Hz, [15](#)
 ADC_REF, [15](#)
 ADC_SAMPLINGRATE, [15](#)
 AttysCommBase, [7](#)
 CHANNEL_DESCRIPTION, [15](#)
 CHANNEL_SHORT_DESCRIPTION, [15](#)

- CHANNEL_UNITS, 16
- closeSocket, 8
- connect, 8
- enableCurrents, 8
- getADCFullScaleRange, 9
- getAccelFullScaleRange, 8
- getAdc_samplingrate_index, 8
- getBiasCurrent, 9
- getBluetoothAdressString, 9
- getBluetoothBinaryAdress, 9
- getIsCharging, 9
- getMagFullScaleRange, 9
- getSampleFromBuffer, 9
- getSamplingRateInHz, 10
- hasActiveConnection, 10
- hasSampleAvailable, 10
- INDEX_Acceleration_X, 16
- INDEX_Acceleration_Y, 16
- INDEX_Acceleration_Z, 16
- INDEX_Analogue_channel_1, 16
- INDEX_Analogue_channel_2, 17
- INDEX_GPIO0, 17
- INDEX_GPIO1, 17
- INDEX_Magnetic_field_X, 17
- INDEX_Magnetic_field_Y, 17
- INDEX_Magnetic_field_Z, 17
- MAG_FULL_SCALE, 17
- MESSAGE_CONNECTED, 17
- MESSAGE_ERROR, 17
- MESSAGE_RECEIVING_DATA, 18
- MESSAGE_RECONNECTED, 18
- MESSAGE_TIMEOUT, 18
- NCHANNELS, 18
- nMem, 18
- oneG, 18
- phys2temperature, 10
- quit, 10
- registerCallback, 10
- registerMessageCallback, 10
- resetRingbuffer, 10
- setAccel_full_scale_index, 11
- setAdc0_gain_index, 11
- setAdc1_gain_index, 11
- setAdc_samplingrate_index, 11
- setBiasCurrent, 11
- start, 11
- unregisterCallback, 11
- unregisterMessageCallback, 12
- AttysCommListener, 18
 - hasSample, 19
- AttysCommMessage, 19
 - hasMessage, 19
- AttysScan, 20
 - getAttysComm, 20
 - getAttysName, 20
 - getNAttysDevices, 21
 - registerCallback, 21
- SCAN_CONNECTERR, 21
- SCAN_CONNECTED, 21
- SCAN_CONNECTING, 22
- SCAN_NODEV, 22
- SCAN_SEARCHING, 22
- SCAN_SOCKETERR, 22
- scan, 21
 - unregisterCallback, 21
- AttysScanListener, 22
- CHANNEL_DESCRIPTION
 - AttysCommBase, 15
- CHANNEL_SHORT_DESCRIPTION
 - AttysCommBase, 15
- CHANNEL_UNITS
 - AttysCommBase, 16
- closeSocket
 - AttysComm, 3
 - AttysCommBase, 8
- connect
 - AttysComm, 3
 - AttysCommBase, 8
- enableCurrents
 - AttysCommBase, 8
- getADCFullScaleRange
 - AttysCommBase, 9
- getAccelFullScaleRange
 - AttysCommBase, 8
- getAdc_samplingrate_index
 - AttysCommBase, 8
- getAttysComm
 - AttysScan, 20
- getAttysName
 - AttysScan, 20
- getBiasCurrent
 - AttysCommBase, 9
- getBluetoothAdressString
 - AttysComm, 4
 - AttysCommBase, 9
- getBluetoothBinaryAdress
 - AttysComm, 4
 - AttysCommBase, 9
- getIsCharging
 - AttysCommBase, 9
- getMagFullScaleRange
 - AttysCommBase, 9
- getNAttysDevices
 - AttysScan, 21
- getSampleFromBuffer
 - AttysCommBase, 9
- getSamplingRateInHz
 - AttysCommBase, 10
- hasActiveConnection
 - AttysCommBase, 10
- hasMessage
 - AttysCommMessage, 19

hasSample
 AttysCommListener, 19
hasSampleAvailable
 AttysCommBase, 10

INDEX_Acceleration_X
 AttysCommBase, 16
INDEX_Acceleration_Y
 AttysCommBase, 16
INDEX_Acceleration_Z
 AttysCommBase, 16
INDEX_Analogue_channel_1
 AttysCommBase, 16
INDEX_Analogue_channel_2
 AttysCommBase, 17
INDEX_GPIO0
 AttysCommBase, 17
INDEX_GPIO1
 AttysCommBase, 17
INDEX_Magnetic_field_X
 AttysCommBase, 17
INDEX_Magnetic_field_Y
 AttysCommBase, 17
INDEX_Magnetic_field_Z
 AttysCommBase, 17

MAG_FULL_SCALE
 AttysCommBase, 17
MAX_ATTYS_DEVS
 AttysScan, 21
MESSAGE_CONNECTED
 AttysCommBase, 17
MESSAGE_ERROR
 AttysCommBase, 17
MESSAGE_RECEIVING_DATA
 AttysCommBase, 18
MESSAGE_RECONNECTED
 AttysCommBase, 18
MESSAGE_TIMEOUT
 AttysCommBase, 18

NCHANNELS
 AttysCommBase, 18
nMem
 AttysCommBase, 18

oneG
 AttysCommBase, 18

phys2temperature
 AttysCommBase, 10

quit
 AttysComm, 4
 AttysCommBase, 10

receptionTimeout
 AttysComm, 4
registerCallback
 AttysCommBase, 10

 AttysScan, 21
registerMessageCallback
 AttysCommBase, 10
resetRingbuffer
 AttysCommBase, 10
run
 AttysComm, 4

SCAN_CONNECTERR
 AttysScan, 21
SCAN_CONNECTED
 AttysScan, 21
SCAN_CONNECTING
 AttysScan, 22
SCAN_NODEV
 AttysScan, 22
SCAN_SEARCHING
 AttysScan, 22
SCAN_SOCKETERR
 AttysScan, 22
scan
 AttysScan, 21
sendInit
 AttysComm, 4
sendSyncCommand
 AttysComm, 5
setAccel_full_scale_index
 AttysCommBase, 11
setAdc0_gain_index
 AttysCommBase, 11
setAdc1_gain_index
 AttysCommBase, 11
setAdc_samplingrate_index
 AttysCommBase, 11
setBiasCurrent
 AttysCommBase, 11
start
 AttysComm, 5
 AttysCommBase, 11

unregisterCallback
 AttysCommBase, 11
 AttysScan, 21
unregisterMessageCallback
 AttysCommBase, 12